Clinitek Status® Analyzer Specifications





Intended Use

The Analyzer is for in vitro use in the semi-quantitative detection of albumin, bilirubin, blood (occult), creatinine, glucose, ketone (acetoacetic acid), leukocytes, nitrite, pH, protein, specific gravity and urobilinogen in urine samples, when Siemens reagent strips for urinalysis are used.

The semi-quantitative calculation of albumin-to-creatinine and protein-to-creatinine ratios in urine samples, when Clinitek® Microalbumin and Multistix PRO® Reagent Strips for Urinalysis are used.

The detection of human Chorionic Gonadotropin (hCG) in urine samples, when Clinitest® hCG test cassettes are used.

Instrument Memory

200 Patient test results, including Operator Identification.

200 Patient details, including name and/or Patient Identification.

Interface	
User:	Touch screen display;
	alphanumeric input; password
	protection.
Computer:	Results are sent to a personal

computer via a RS-232 serial

port using a 9-pin null modem

Power Requirements

- 110V AC +/- 20% 45-65 Hz.
- Battery powered operation (optional) 6 AA non-rechargeable alkaline batteries, not included.

Dimensions/Weight

Depth:	10.7 inches Width: 6.7 inches
Height:	6.2 inches Weight: 3.65 lb.
	Clinitek Status Analyzer only
	without batteries or power supply.

Operation

Temperature Range:	41°F - 104°F
Humidity Range:	20% - 80% Relative Humidity,
	non-condensing.
Calibration:	Automatic, self-calibrating.

Urinalysis has become simpler and more reliable with the Clinitek Status Analyzer.

Answers for life.



Clinitek Status Analyzer

Ordering Information

Uninghasia Chuina	Product Code	Packaging Configurations	Glucose	Bilirubin	Ketone	Specific Gravity	Blood	Hd	Protein	Urobilinogen	Nitrite	Leukocyte	Albumin	Creatinine	hCG
Urinalysis Strips															
Multistix® 10 SG	2161	100	•	•	•	•	•	•	•	•	•	•			
Multistix® 8 SG	2164	100	•		•			•	•		•				
Multistix® 7	2165	100	•		•		•	•	•		•	•			
Multistix® 5 Multistix® SG	2389	100							•		•	•			
N-Multistix® SG	2177	100	•												
Uristix®	2176	100		_	_			_		_					
Uristix® 4	2166	100	•												
Hema-Combistix®	2182	100	•				•	•	•						
Urine Reagent Stri			tic	,											
Multistix PRO® 10LS	1554	100	•						•		•	•		•	
Multistix PRO 10LS	1566	25	•		•	•	•	•	•		•	•		•	
Clinitek® Microalbumin	2083	25											•	•	
Cassette		l													
Clinitest® hCG	1760	25													٠
Quality Controls	1700	23													•
Chek-Stix® Urinalysis															
Control Strips	1360	25													
(Positive)	1300	23					Ĭ								
Chek-Stix® Combo Pak	1364	25													
(Positive/Negative)			•	•	•	•	•	•	•	•	•	•			
Analyzer															
Clinitek Status®	1735	Each													
Climitals Duintan Danas															-
Clinitek Printer Paper	5773	5													

Table of Results ▶

The results shown in shaded areas will be marked as positives, if "mark positive results" is selected in Instrument Set Up.
The results will be marked by asterisks when displayed, when printed and when the data is transferred to a host computer.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Table of Results

		Reager	nt Strip T	ests				
Test	Abbreviation	Units		Reported	Results			
rest	Appreviation	Units	Normal	System	Plus S	ystem		
			Negative	500	Negative	2+		
Glucose	GLU	mg/dL	100	>=1000	Trace	3+		
			250		1+			
Bilirubin	BIL		Negative	Moderate	Negative	2+		
Dilliabili	DIL		Small	Large	1+	3+		
			Negative	40	Negative	2+		
Ketone	KET	mg/dL	Trace	80	Trace	3+		
			15	>=160	1+	4+		
Specific			<=1.005	1.020				
Gravity	SG		1.010	1.025	No diff	erence		
Gravity			1.015	>=1.030				
Occult			Negative	Small	Negative	1+		
Blood	BLO		Trace-lysed	Moderate	Trace-lysed	2+		
Diood			Trace-intact	Large	Trace-intact	3+		
			5.0 6.5	8.0				
рН	рН		5.5 7.0	8.5	No difference			
			6.0 7.5	>=9.0				
Protein	PRO		Negative	100	Negative	2+		
Multistix PRO		mg/dL	15	300	Low	3+		
			30		1+			
Protein	PRO	mg/dL	Negative	100	Negative	2+ 3+		
All other urinalysis			Trace	>=300	>=300 Trace			
strips			30		1+			
		E.U./dL	0.2	4.0				
Urobilinogen	URO		1.0	>=8.0	No diff	erence		
			2.0					
Nitrite	NIT		Negative	Positive	No diff			
			Negative	Moderate	Negative	2+		
Leukocytes	LEU		Trace	Large	Trace	3+		
			Small		1+			
Albumin	ALB	mg/L	10	80	No diff	erence		
			30	150				
c	CRE		10	200				
Creatinine		mg/dL	50	300	No diff	erence		
			100	. 200 11:				
Albumin:	A:C	,	<30 Normal	>300 High Abnormal				
Creatinine*		mg/g	30 - 300		No diff	erence		
			Abnormal					
Protein:			Normal	300				
Creatinine	P:C	mg/g	Dilute Normal	Abnormal	No diff	erence		
		iligig	150	`E00	No difference			
			Abnormal	>500 Abnormal				
		Cass	ette Tes					
			hCG					
Ulara			Negative	hCG Positive				
Human Chorionic	hCG	N/A	hCG		No diff	difference		
Gonadotropin	iica	11// (Borderline		No difference			
			Repeat in 48-72 hours					

^{*}Positive Albumin: Creatinine ratios are not marked when using Clinitek Microalbumin Reagent Strips for Urinalysis.

Siemens Healthcare Diagnostics Inc. 1717 Deerfield Road Deerfield, IL 60015-0778 USA

© 2008, Siemens Healthcare Diagnostics Inc. Order No. A91DX-0351061-C2-4A00 Printed in USA

Clinitek Status, Clinitek, Clinitest, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners.