



Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

PROFESSIONAL MEDICAL PRODUCTS

PC-900B CAPNOGRAPH AND OXIMETER

User Manual

REF PC-900B (Gima 33698)



Shenzhen Creative Industry Co., Ltd.
Floor 5, BLD 9, BaiWangxin High-Tech Industrial Park,
Songbai Road, Xili Street, Nanshan District,
518110 Shenzhen, P.R. China
Made in China



Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany



Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



Notice

Thanks for buying PC -900B Capnograph and Oximeter.

This manual is copyright reserved. It is prohibited to copy, duplicate or translate into other languages without our written permission.

Please read this manual carefully and then follow its instructions when operating this monitor.

It is not permitted to open the monitor's main cover, modify or disassemble it without our permission or official service training.

The buyer will not be advised of technology updates which do not influence the monitor's key functionality. Furthermore, please pay attention to the difference between the parts or components provided as information in this manual.

Shenzhen Creative Industry Co., Ltd.

Floor 5, BLD 9, Baiwangxin High-Tech Industrial Park,

Songbai Road, Xili Street, Nanshan District,

518110, Shenzhen, China

CONTENTS

1 Preface	1
1.1 Brief	1
1.2 Warranty and Maintenance	1
1.3 Safety Requirements	2
2 Technical specifications and characteristics	3
3 Introduction of Monitor	6
4 Patient connection	8
4.1 CO ₂ and Respiration rate measurement	8
4.2 CO ₂ Measure principle	10
4.3 Oximeter density measurement (optional)	11
4.4 Notice	12
5 Screen display and Operation	15
5.1. Screen main display menu	15
5.2 Initial Monitoring Screen	17
5.3 The Main Menu	18
5.4 CO ₂ SET Menu	19
5.5 SpO ₂ SET Menu	21
5.6. TIME SET Menu	22
5.7. Sound SET Menu	24
5.8. Trend	24
5.9. NEW PATIENT Menu	27
6 Charging, Maintenance, Cleaning	29
6.1 Charging	29
6.2 Maintenance	29
6.3 Cleaning	30
7 Trouble Shooting Analysis	32
Appendix 1. Explanations of Terms in this Manual	34
Appendix 2. Changing compensation of balance gas	35
Appendix 3. Calibration of EtCO ₂ Accuracy	37
Appendix 4. Part Numbers and Consumables listing	41
Appendix 5. Guidance and manufacturer's declaration -Electromagnetic compatibility	43

1 Preface

1.1 Brief

The purpose of this manual is to provide the user with a brief understanding of the characteristics, functions and operation of the monitor thereby preventing incorrect operation and user error.

This monitor can measure four physical patient/animal parameters at the same time: concentration of EtCO₂, respiration rate, heart pulse rate and saturation of SpO₂ (optional). The monitor you bought may have two or more functions mentioned above but this manual can be used in common for the applicable functions.

1.2 Warranty and Maintenance

Warranty

This monitor has a warranty of 12 months from the date of purchase. Reusable SpO₂ sensors and the battery included have a 12 month warranty. All other accessories have a warranty of 3 months or an “out of box” warranty for disposable items.

The following will invalidate the warranty:

- if the monitor is damaged due to misuse or incorrect operation (i.e. without following the user manual instruction)
- the monitor is damaged due to incorrect connection with another instrument
- the monitor is accidentally damaged or dropped
- if the user modifies or changes the monitor without written authority of the company
- if the serial number is deliberately damaged, torn off or unreadable.

Maintenance

If the monitor is non-functional outside of the warranty period, the manufacturer or distributor will offer an estimate for repair. The maintenance, repair or calibration would be carried out at local distributor, unless detailed in a

specific written agreement.

Re-packing for Repair or Calibration

It is recommended to use the original packing boxes and packing materials when returning for repair or maintenance.

1.3 Safety Requirements

For the purposes of safety, please read the following and abide by these instructions for medical instrumental products.

Warning: Indicating the possible injury on patient or operator.

- This monitor is not MRI compatible and is not suitable for use within the magnetic field during the operation of MRI or CT. However, the sample lines supplied alongside the unit by the distributor are MRI compatible and may be extended into the MR or CT field. In this case, the monitor must remain outside of the room.
- The use of accessories and cable other than those specified, with the exception of cables sold by the manufacturer of the device as replacement parts for internal component, may result in increased emissions or decreased accuracy of the device.
- Only use manufacturer designated accessories to ensure compliance with appropriate standards
- It is not allowed to remove the cover of the monitor.
- This monitor provides concentration of EtCO₂, respiration rate, oxygen saturation and pulse rate. This data only provides assistance for diagnosis and actual diagnosis shall be made by suitably qualified clinical staff using all the clinical information and symptoms.
- In order to prevent pressure sores and correct circulation the SpO₂ sensor must be repositioned regularly, depending on the type of sensor used.

2 Technical specifications and characteristics

Intended Use

The Capnograph and Oximeter is designed for monitoring the vital physiological signs of the patient. It is used for non-invasive continuous monitoring of oxygen saturation (SpO₂), pulse rate, CO₂ and respiration rate.

The Capnograph and Oximeter is intended for use in adults, pediatrics and infants in a hospital environment. It is intended to be used only under regular supervision of clinical personnel.

EtCO₂

Sampling Mode: Sidestream

Method: Non-dispersive Infrared Spectroscopy

Range: 0 – 150mmHg or 0 – 20kPa or 0 – 20% (v/v)

Accuracy: ±2mmHg for EtCO₂ range 0 - 40mmHg

±5% for EtCO₂ range from 41 - 70mmHg

±8% for EtCO₂ range from 71 - 100mmHg

±10% Over 100mmHg

Note: The accuracy of CO₂ concentration measurement is influenced by any interfering gas and/or vapour, for example N₂O gas can raise the CO₂ reading (2-10%), and Helium and O₂ can reduce the CO₂ reading (1-10%), so compensation should be set in the balance gas MENU to meet the accuracy requirements if such gases or vapours are present.

Update/Averaging Time: Option of every breath or 10, 20 or 30 seconds

Warm Up Time: <20 seconds

Sample Flow Rate: 50–250ml/min User Adjustable. Default=100ml/min

Memory: 24 hours on Screen Trend and Numeric

Sensor: <25g Single Use Gas Sample Line and Adaptor for Intubated and /or Non Intubated Patients

Respiration Rate

Range: 3 - 150 breaths/minute

Accuracy: ±1% of reading or ±1 breaths/min whichever is greater

Memory: 24 hours on Screen Trend and Numeric

SpO₂ (optional)

Transducer: Dual-wavelength LED

Range: 0 - 100%

Accuracy: ±2% for SpO₂ range from 70 - 100%,
±3% for SpO₂ range from 50 - 69%

Memory: 24 hours on Screen Trend and Numeric

Patient Modes: Adult and Pediatric

Pulse Rate (optional)

Range: 30 – 250bpm

Accuracy: ±2% for PR range from 30 - 250bpm

Memory: 24 hours on Screen Trend and Numeric

Power

AC Input: 100V - 250V, 50Hz/60 Hz to 5VDC Adapter with
5V mini USB adapter Cable.
Optional Vehicle 12V to 5V Mini USB Charger Lead.

Battery

Type: Built-in rechargeable lithium battery pack ,(3.6V, 3000mAH)

Charging Time: 4 hours from flat

Operating Time: 10 hours on full charge

Operating Conditions

Temperature: +5 to +50°C

Humidity: < 93% % (non-condensing) = < 29.45 hPa

Atmospheric pressure: 70 - 120 kPa

Storage Conditions

Temperature: -30to +70°C

Relative Humidity: <93% (non-condensing)

Atmospheric pressure: 50 - 120 kPa

Dimensions of Monitor

Size: 70 x 160 x 40mm (W x H x D)

Weight: Monitor 380g, Weight on Airway ETT/LMA <25g.

Warranty & Maintenance/ Calibration

One year warranty on main unit and lithium ion rechargeable battery

Auto self-zeroing calibration, annual calibration check recommended

IP rating

IP32 when used in specified carry case.

CE & Product classification

As per IEC 60601- 1 / CSA601.1 / UL2601-1

Type of Protection

Class II (When used with UK/EU Power Supplies)

Degree of Protection: Type BF-Applied Part

Mode of Operation: Continuous

Electro-Magnetic Compatibility: Group I, Class A

93/42/EEC Medical Device Directive Compliant

EC-Representative:

Shanghai International Holding Corp. GmbH (Europe)






Eiffestraße 80, 20537 Hamburg Germany

CE 0123

3 Introduction of Monitor



Figure 3.1

- (1) Screen: Displays waves, menu, alarm and all measuring parameters.
- (2)  / : Function button:
 - ▲ a) When menu (except the TREND menu) is activated, press this button to move the cursor.
 - b) When the TREND menu is activated, this button changes between the trend graph and data table
- (3) : On the main display, to press this button to silence alarms for 2 minutes.
- (3) : Press this button to move the cursor when menu is activated.
- (4) : Multifunction button.
 - a) Press this button to increase figures on the menu.

b) In the main display screen, press this button to freeze the display waveform (if frozen, the data which prints will be that shown on the screen).


(5) -: Press this button to decrease figures.

(6) **ENTER**: Confirmation button;

a) Press this button to "Confirm" on the menu.

b) In the main menu, press this button restart the pump if it has automatically switched OFF.

c) If the device is connecting with Bluetooth printer, press this button for 2 seconds to print capnography and other result parameters (EtCO₂, RR, SpO₂, PR).

(7)  Press this button to enter or quit menu or change display

(8)  Power button: hold for >2 seconds to activate

(9) Indicator **POWER**: Blue LED is lit when the Monitor is either switched ON or subject to external power when not switched ON.

If the yellow LED is lit, the internal battery is being charged.

(10) **CO₂**: The faucet of filter, blue color indicator flashes if the filter is off. When the filter is plugged in, the indicator color will change to blue, and it will change to red during occlusion or pump err.

(11) **SpO₂**: The socket of SpO₂ (optional).

(12) DC5V Mini USB Charging interface. Note: this interface must only be connected to a device which meets safety standards.

(13) Exhaust outlet: Do not occlude.

(14) Speaker location

(15) Battery Compartment with clip on Battery Door

(16) Hanging Point for Lanyard if required.

4 Patient connection

4.1 CO₂ and Respiration rate measurement

Push in and twist 45° clockwise to connect the Filter/Water Trap to the Connector on the top of the Monitor. Attach the selected Gas Sampling Line to the CO₂ filter/Water Trap Female Luer Connector (Use a Male to Male Luer adapter if necessary) and then select a sampling point as close as possible to either the Animal or the Ventilator Breathing Circuit.




Figure 4.1

⚠ WARNING ⚠ :

Do not use the Monitor if the filter/water trap is not installed to avoid contamination and damage to the IR measurement cell.

In order to avoid vapor and respiratory mucus entering into the IR Cell, the machine must be used with the Filter/Water Trap .




Figure 4.2

Instruction for use of the Filter/Water Trap

- 1.) Insert the convex cleat of Filter/water Trap into the notch of the inlet port of device and turn 45° clockwise.
- 2.) Attach male luer lock sample line connector to Filter/Water Trap (Use a Male to Male Luer adapter if the sample line has Female Luer connector)
- 3.) Connect the other end of the Sample line to the chosen sampling point of Animal or Ventilator Circuit.
- 4.) Change the Filter/Water Trap as needed. If the Filter/Water Trap becomes dirty or the occlusion alarm is activated when it is dry then the Filter/Water Trap must be replaced.

 WARNING  :

Ensure that connections are air tight as if there is leakage, measured values are likely to be inaccurate.

 WARNING  :

Use only recommended original bespoke Filter Water Trap to ensure accuracy.

4.2 CO2 Measure principle

1. Theory introduction

The device working theory is NON-DISPERSIVE INFRA GAS ANALYZER. The device has an AUTO ZERO ADJUSTMENT SYSTEM and GAIN CONTROL

The principle is based on the fact that CO₂ molecules absorb infrared light energy of specific wavelengths, with the amount of energy absorbed being directly related to the CO₂ concentration. When an IR light beam is passed through a gas sample containing CO₂, the electronic signal from an **infrared sensor**(which measures the remaining light energy), can be obtained. This signal is then compared to the energy of the IR source, and calibrated to accurately reflect CO₂ concentration in the sample. To calibrated, the infrared sensor's response to a known concentration of CO₂ is stored in the monitor's memory.

Further, on the passage of the sample gas with a three-way valve, with the change of temperature and time, the valve leading to the pure air will take 3-4 seconds in order to adjust the zero point.

In addition, the circuit module has the **atmospheres absolute pressure sensors** and the **flow measurement and control of pressure sensor**. Modules can measure atmospheric pressure, and atmospheric can compensate the calculation for the concentrations of carbon dioxide which improve the design accuracy.

Then the monitor (CO₂ module) determines CO₂ concentration in the breathing gases by measuring the amount of light absorbed by these gases. EtCO₂ is display as a numerical value in millimeters of mercury (mmHg), percent (%), or kilopascals (kPa). In addition, a CO₂ waveform (capnogram) may be displayed which is a valuable clinical tool that can be used to assess patient airway integrity and proper endotracheal tube placement. Respiration rate is calculated by measuring the time interval between detected breaths.

2. Automatic Offset Calibrations:

The device was designed to automatically perform calibrations in order to correct for changes in temperature, altitude and electronic component drift. The air surrounding the device may have elevated concentrations of CO₂ present (such as in an enclosed compartment or room with poor ventilation). Therefore, we recommend use in well ventilated locations to ensure that the CO₂ baseline does not cause inaccuracy.

3. The Moisture Separation System:

This instrument uses a patented filter/water trap which can filter a large amount of moisture whilst maintaining a minimum dead space thereby improving the accuracy of the waveform. Please note that if the Filter/Water Trap becomes full of water or dirt the display will show "OCCLUSION", the operator needs to change the filter/water trap. The old filter can be reused after natural drying in a ventilated and dry environment. Discard the old filter/water trap if dirty.

4.3 Oximeter density measurement (optional)

Theory introduction

SpO₂ is measured by Pulsating oximetry. This is a continuous, non-invasive method to measuring hemoglobin oxygenation saturation. It is determined the number of sensor light emitted from the light source side penetrating the animal tissue (such as a finger or ear), to the receiver sensor. Sensors measure the wavelength of the red LED is typically 660nm, infrared LED is 940nm. The maximum output power of the optional LED is 4mW.

The amount of light passing through depends on several factors, most of which is constant. However, the arterial blood flow that is one of these factors varies with time, because it is pulsating. By measuring the light absorption of the pulsation period, it is possible to obtain the oxygen saturation of arterial blood. Detect pulsating itself can give a "plethysmography" wave and pulse rate signal.

It is also recommended to use Pulse Oximetry for ventilated or sedated animals. Measurement will begin when a finger is put into the sensor clip, meanwhile, the photoplethysmogram wave will appear on the screen, after several seconds

the oxygen saturation and pulse rate appear. The monitor will give a pulse tone sound when each heart beat happens. The tone will change to an alarm tone if the values of SpO₂ and Pulse Rate breach the alarm level settings. The volume of pulse beep can be adjusted by the item **BEEP VOLUME** in the SOUND SET menu. The pulse beep tone will disappear under the silent condition.




Figure 4.3

The use of different SpO₂ sensors

There are a number of different SpO₂ Sensors for use with this monitor. Please see brochure or listing at rear of Manual for details.

PLEASE NOTE: when SpO₂ is not being monitored the probe should be disconnected from the monitor to save battery life, or the two windows of sensor should be kept face to face, otherwise the light window will remain operational and the photoplethysmogram wave will be disordered and the screen will display “**FAIL SEARCH**”.

4.4 Notice

1. Caution:

Conditions of electromagnetic influence, for example: electrosurgical devices, MRI, CT etc., may cause incorrect operation.

This device is not MRI/CT Compatible.

The filter/water trap should be taken off and replaced when it is nearly full of water, otherwise water ingress may cause irreversible damage for IR measurement

detector cell. Be sure that the collecting pipe is not occluded to avoid stressing the inner sampling pump and reduction of pump life.

2. Attention: other important information.

1.) CO₂:

The approved sampling lines provided by or specified by the manufacturer or distributor, shall be used, otherwise readings may be inaccurate.

Fast changes in ambient Temperature may cause inaccuracy and in this instance the Display will show "TEMP IMBALANCE".

The measured data may be influenced by different kinds of anaesthetic gases. If it is required to calibrate interference gases please refer to Appendix 2.

Any circumstances of blocking of the gas sampling line, such as bending, folding, contamination blocking the sampling tube and filter or water trap etc. may lead to inaccurate measurement.

Serious respiratory conditions leading to exhaled CO₂ concentration being extremely high or low, e.g. EtCO₂ lower than 0.5% or higher than 11%, may generate inaccurate measurement.

Any air leaks in the sampling line circuit will seriously influence accuracy of data measured and waveform shape.

2.) Oximeter:

The monitor's measurement of SpO₂ may be influenced by strong ambient light. Therefore the user should unplug the SpO₂ Sensor when it is not being used.

Accuracy of oximeter readings will be influenced if there is imaging dye in the blood or if CO has been inhaled by the Patient.

Only use original SpO₂ probes approved for use with this Monitor.

Always make sure that the sensor is not contaminated or broken before use.

Always take care to check that the sensor is applied correctly.

Warning:

Do not use the SpO₂ sensor if it is damaged or dirty.

If shock, low blood pressure, serious blood vessel constriction, serious anemia,

very body low temperature, artery block near sensor or incomplete heart asystole occur the pulse signal may disappear.

5 Screen display and Operation

5.1. Screen main display menu








Figure 5.1

1. The first line of data shows time (hour, minute)/patient ID, animal type: THICK or THIN tissue, the memory area full indicator (☐), alarm sound closing (!🔔), silence (🔕) or non silence(🔔), Bluetooth symbol (📶) and battery indicator 🔋.


Attention:

- When the memory full indicator is displayed, further data cannot be stored. If you want to save the new data effectively, you need to enter the NEW

PATIENT menu to delete the data in the storage area, or to change patient ID. Alternatively, select AUTO LOOP to overwrite the oldest data when memory is full, please see the details in 5.9 NEW PATIENT

b) If the symbol  appears, the menu is locked, the setting menu will be disabled unless user press the three buttons  ,  ,  -at the same time, or


enters engineer menu to unlock the menu (Refer to Appendix 2. ENGINEER MENU: Changing compensation of balance gas)

c) The middle part of the screen shows results data: EtCO₂ concentration, respiratory rate, inhaling CO₂ concentration (optional), oxygen PLETH, exhaling or inhaling state (during exhaling,  becomes blue color).

The bottom area shows CO₂ respiratory wave. If it is equipped with SpO₂, it will show SpO₂, pulse, oxygen PLETH waveform and histogram. When the pump is not operating "PUMP OFF" will appear on the screen. If the filter/water trap is NOT inserted into the inlet port, the screen will show 'LINE OFF, the pump will also be automatically switched off to prevent ingress to the unprotected IR detector cell.

Alarm indication:

1.) If the EtCO₂'s value exceeds the limit of high or low alarm level, the word 'EtCO₂' will flash and alert with the audible high priority alarm. This high priority alarm will also sound for respiration rate, SpO₂ and pulse rate alarms.

2.) If the battery level is almost fully depleted the battery  indicates completely empty, the monitor will alarm continuously and will shut down automatically.


3.) When the apnea alarm is turned on and apnea occurs the monitor will give a high priority audio/visual alarm. The screen will flash the message 'APNEA' (meaning no EtCO₂ has been detected for a certain time period) and if the sound alarm is turned on, it will alert a high priority audible alarm.

4.) When the SpO₂ sensor is disconnected or not applied, the screen will flash the message 'SENSOR OFF'. If a heart beat pulse is not detected for a period of time, the screen will flash the message 'FAIL SEARCH'.

5.) The volume of continuous or interval alarm tone sounds mentioned above can be adjusted up and down by the menu item **ALARM_VOLUME**. The sound will be inaudible under the silent condition. If the alarm volume is 0, the silence indicator in the main menu will show '!'.

6.) All the parameter alarms for over limits and apnea alarm, will lead to the flashing of the red alarm indicator on the panel.

5.2 Initial Monitoring Screen

Long press (about 3 seconds) power key "  " to start the monitor, the initial monitoring screen is as shown below:




Figure 5.2

In this menu, press **▲ / +** button or **▼ / -** button to move the cursor, then press the **ENTER** button to select YES or NO. If selecting "YES" then the monitor enters the New Patient menu directly. If selecting "NO" or there is no

any operation in 8 seconds, then the monitor enters main display screen. To disable this prompt, enter the New Patient menu screen. If "POWER ON ID PROMPT" is set as "NO", the monitor will disregard the initial monitoring screen (see figure 5.2) and enter into main display screen directly (refer to Section 5.9 NEW PATIENT MENU for details).

5.3 The Main Menu







Figure 5.3

Press the MENU button  to enter the Main Menu to set monitor parameters

 **WARNING** : **All Menu Settings are LATCHING and remain when the Monitor is powered off. Ensure that all necessary settings are reviewed and are suitable for the patient BEFORE use.**

This menu includes the following options:

The setting menu for CO₂: **CO2_SETUP**

The setting menu for SpO₂: **SPO2_SETUP**

The trend menu: **TREND**

The time menu: **TIME_SETUP**

The sound menu: **SOUND_SETUP**

The new patient menu: **NEW PATIENT**.

In this menu, to press ▲ or ▼ button to move the cursor up or down to highlight an option and Press the ENTER button to select and enter the next level of the menu. To return to the Main menu select EXIT option and press ENTER (not available on Trend screen).

5.4 CO₂ SET Menu

CO2 SET		
EiCO2	ALARM_H	50.0mmHg
	ALARM_L	19.0
RESP	ALARM_H	30RPM
	ALARM_L	08
FLOW RATE	100	ml/Min
APNEA TIME	30	S
UNIT		mmHg
CO2 PUMP		ON
AUTO OFF TIME	10	Min
SWEEP SPEED		FAST
WAVE SCALE	54	mmHg
EiCO2 AVER		1BREATH
LOAD DEFAULTS		
EXIT		

Figure 5.4

In this menu, press ▲ or ▼ button to move the cursor up or down, press +

button or - button to change the data highlighted by the cursor.

To return to the main menu highlight EXIT and press the **ENTER** button. If you want to return the monitor to its default settings highlight LOAD DEFAULTS and Press the **ENTER** button.

This menu includes the following setups:

- 1.)The high alarm limits of EtCO₂:**EtCO2 ALARM_H**: 22-99mmHg, off
- 2.)The low alarm limits of EtCO₂:**EtCO2 ALARM_L**: off, 10-60mmHg
- 3.)The high alarm limits of respiration rate: **RESP ALARM_H**:5-60 breaths/min, off
- 4.)The low alarm limits of respiration rate: **RESP ALARM_L**: off, 4-40 breaths/min
- 5.) Pump flow rate setup: **FLOW-SET**: 50 -250ml/min
- 6.)The setup of apnea time: **APNEA TIME**: 15s-44s, off
- 7.) The unit of CO₂:**CO2 UNIT** : %, mmHg or kPa
8.)Pump switch: ON or Off
- 9.) Pump auto-closing time: **AUTO-OFF-TIME**: 10-30min
- 10.)Screen speed of capnograph : **SWEEP SPEED**: SLOW, NORMAL or FAST
- 11.) CO₂ Wave scale: **WAVE SCALE**: 54mmHG or 76mmHG
- 12.) EtCO₂ average computation time: **EtCO₂ Averaging**: every breath, 10sec, 20sec, 30sec
- 13.) Default reload: **LOAD-DEFAULTS**
- 14.) Exit: **EXIT**

Attention:

Pump auto-closing time means that the pump will automatically be closed down when no respiration occurs in the set period (default 10 min).

The wave scale means the maximum value of waveform amplitude display but it does not mean data on full-scale. Data on full-scale still means 99mmHg.

Default values are as follows:

EtCO₂ alarm high limit: 50 mmHg
 EtCO₂ alarm low limit: 19 mmHg
 RESP alarm high limit: 30 breaths/min
 RESP alarm low limit: 08 breaths/min
 FiCO₂ alarm high limit: OFF
 FLOW_SET: 100 CC/Min
 Apnea time: 30S
 CO₂ unit: %
 CO₂_PUMP: ON

AUTO_OFF_TIME: 10 Min
SWEEP SPEED: NORMAL
EtCO₂ Averaging: 1 Breath
WAVE SCALE: 54mmHg

5.5 SpO₂ SET Menu




Figure 5.5

In this menu, press ▲ or ▼ button to move the cursor up or down, press + button or - button to change the data highlighted by the cursor.

To return to the main menu highlight EXIT and press the **ENTER** button.
If you want to return the monitor to its default settings highlight LOAD DEFAULTS and Press the **ENTER** button.

This menu includes the following setups:

- 1.)The low alarm limits of SpO₂: **SPO2 ALARM_L**: off, 50%-99%
- 2.)The high alarm limits of pulse rate:
P_RATE ALARM_H: **50-350**beats/m, OFF
- 3.)The low alarm limits of pulse rate:
P_RATE ALARM_L: OFF, **20-150**beats/m

4.)Wave curve selection: **CURVEWAVE**: FILL or LINE

5.)Renewing of defaults. **LOAD DEFAULTS**

The wave curve selection means that: FILL indicates the beneath part of photoplethsmogram is filled. LINE indicates the photoplethsmogram is drawn in curve line.

Default values as follow:

SpO₂ alarm low limit: 92%

Pulse Rate alarm high limit: 130bpm

Pulse Rate alarm low limit: 50bpm

Curve: Line

5.6. TIME SET Menu

TIME SET	
YEAR	13
MONTH	01
DATE	10
HOUR	21
MINUTE	18
SAVE	
EXIT	

Figure 5.6

In this menu, press ▲ or ▼ button to move the cursor up or down, press +

button or - button to change the data highlighted by the cursor.

Attention:

Any time adjustment will delete any stored trend data, so please take care before making this adjustment.

The procedure is as follows:

- 1.) Change time.
- 2.) Move the cursor to *SAVE* then press the ENTER button to enter the confirm menu, as the below figure showing.
- 3.) YES is already selected (highlighted in white) and if you wish to confirm this change press Enter if you do not wish to confirm the change move the cursor and highlight NO and press Enter.
- 4.) Only by confirming can the time adjustments be made.




Figure 5.7

5.7. Sound SET Menu




Figure 5.8

In this menu, press ▲ or ▼ button to move the cursor up or down, press +

button or - button to change the data highlighted by the cursor.

This menu includes following setups:

Pulse sound volume: **BEEP_VOLUME**: 0(OFF)-8

Alarm sound volume: **ALARM_VOLUME**: 0(OFF)-8

If the alarm volume is 0, the top line in the main menu will show '!

5.8. Trend

The graph trend





Figure 5.9

The monitor stores EtCO₂, PR, SpO₂ and PR as a group of data every 12seconds (Adjustable in Store Interval under New Patient menu) with accumulated trend up to 24hours respectively. The stored data is retained even the device is shut off.

The symbol  will appear on screen when the storage is full. There are three options to further store the data.

- 1.) Change patient ID under NEW PATIENT menu.
- 2.) Change store mode to AUTO LOOP under NEW PATIENT menu, in auto loop mode new data will be stored and overwrite old data when reaches its limits.
- 3.) Select CLEAR MEMORY under NEW PATIENT menu to empty the stored data.

This figure shows that the time base for the trend page is 1 hours and every point indicates the result of every 12 second. The top line of this page indicates patient's ID number, the start time of this page (date/month/year hour : minute), current page no. and sum pages (24 pages in total).

If in the corresponding time to the one page of trend table, the user turns off and turns on the device once or more times the trend table will show one or several blue vertical lines with full amplitude, at this time press ▼, then the top row will

display the initial information at that turn on time: patient's ID number and initial time. The correspondingly initial blue vertical line will become white one. Press ▼ again, the second initial time will display (if turned off and on for several animals).

The time at beginning and ending parts of abscissa in this picture respectively indicates the beginning and ending time for trend of this page.

If the data is not complete, it shows the monitor was turned off although it has not completed 2 hours' record.

In this menu, press ENTER button to change the trends of CO₂ concentration, respiration rate, SpO₂ and pulse (the latter 2 parameters are selectable).

In this menu, press + button or - button to change the page of trend.

In this menu, press ~~X~~ / ▲ button to change graph trend to table trend .

In this menu, press MENU button to quit this menu and return to the main display.

The table trend

the table No.in one page	patient ID	the beginning time of the large page(one hour):	white color			
time (hour:min:sec)	ID	06/02/2014	07:12	09/23	the page No./the pages sum	the abbreviation of parameter
	00/14	Et	PR	SPO2	PR	
	07:28:12	00	00	00	00	
	07:28:24	00	00	00	00	
	07:28:36	00	00	00	00	
	07:28:48	36	12	99	78	
	07:29:00	38	12	98	70	
	07:29:12	36	12	98	76	
	07:29:24	37	12	98	77	
	07:29:36	39	12	98	77	
	07:29:48	37	12	98	76	
	07:30:00	36	12	98	79	
	07:30:12	36	12	99	78	
	07:30:24	38	12	98	70	
	07:30:36	36	12	98	76	
	07:30:48	37	12	98	77	
	07:31:12	39	12	98	77	
	07:31:24	37	12	98	76	
	07:31:36	36	12	98	79	
	07:31:48	39	12	98	77	
	07:32:00	37	12	98	76	
	07:32:12	36	12	98	79	

the parameter data results: if all zero blue color; otherwise, green color

Figure 5.10

In this graph trend menu, press ~~⏏~~ / ▲ button to change graph trend to table trend. Press ~~⏏~~ / ▲ button again, to return to graph trend .

Every trend table shows **20** groups of data, including time, EtCO₂ (Et), respiration rate (RR), SpO₂, pulse rate (PR). The store interval is adjustable at 12seconds in STORE INTERVAL under NEW PATIENT menu.

There are 24 sum pages when the storage is full. Each page contains **15** trend table and each trend table contains **20** groups data. The **15** trend table in one page can be reviewed by ▼ button. The table no. is indicated on left top of the screen as above figure showing.

In fully stored status, 24 pages can be paged up or down by + button or - button.

The page no. is indicated on right top of the screen as above figure showing.

To quickly check if the four parameters of a data group are all zero, the display will display the parameter columns in blue.

5.9. NEW PATIENT Menu




Figure 5.11

In this menu, press ▲ or ▼ button to move the cursor up or down, press + button or - button to change the data highlighted by the cursor.

Press MENU button, then to exit this menu and enter the main menu.

This menu includes the following setups:

- 1.) **CLEAR MEMORY:** to delete all the historical data so as to store new data
- 2.) **MEM MODE:** to change store mode between manual data deletion (STOP WHEN FULL) and automatic overwriting of the oldest data (AUTO LOOP).
- 3.) **ID: patient's ID,** press "Enter" key to enter or exit from the Set menu. Press

+ button or - button to move the cursor up or down, press ▲ or ▼ button to


change the data highlighted by the cursor.

- 4.) **TYPE:** animal type, THICK or THIN options
- 5.) **STORE INTERVAL:** adjustable at 4/6/12 seconds
- 6.) **POWER ON ID PROMPT:** to set if the monitor enters into the "input new patient" menu when power on the monitor.
- 7.) **SAVE:** store the changes made (it needs to be confirmed by the new menu due to possibly substitution to the original data of the same ID of patient)
- 8.) **EXIT:** to quit the current menu but not to store any changes to the setup

6 Charging, Maintenance, Cleaning

6.1 Charging

Connect the AC/DC power adapter via the Mini USB port turn on the unit. The unit will charge the battery with power at the same time as operating. The battery charge will end after battery is full.

The battery of this unit is a long life rechargeable lithium battery. When the unit is operated on battery only the battery indicator shows the battery's charge level on the screen. When the battery charge level is low, the indicator will flash red , and the external 5VDC power must be connected as soon as possible.

After DC power is connected, the monitor will recharge the battery, and will stop charging after the battery has fully recharged. Operation time for a fully charged unit is > 10 Hours. Charge time is approx. 4 Hours.

Battery replacement method:

Note that the operation must be done with the DC Charger disconnected ensuring that the operator's safety is not compromised.

Press down and slide off to remove the battery cover, then carefully disconnect and remove the battery. Reverse this procedure to replace the new battery and re-fit the battery cover.

NOTE: Any battery that is removed and no longer required must be properly disposed of by following national and local regulations.

6.2 Maintenance

If the monitor appears abnormal (e.g. software system is halted), then to reboot the device hold the Power ON/OFF button down for 5 seconds.

OCCLUSION: If the Display shows 'occlusion', check if the filter /water trap and/or sampling line tubing or connectors are blocked. Replace as necessary and clear the occlusion or switch OFF to prevent damage to the sampling pump.

Please do not let alcohol, cleaning reagent or sterilizing reagent into filter/water trap. Check that the filter/water trap is dry and clean before it is used. Replace the filter/water trap if it is dirty, shows any sign of contamination or if in any doubt about its condition.

6.3 Cleaning

Warning: Before cleaning the monitor and probe, turn off power and remove from any charging source.

1.)Cleaning the Monitor

It is recommended that the Monitor is used in the supplied Carry Case which offers protection from both contamination, liquid ingress and damage.

Do not sterilize by high pressure, autoclave or washer

Do not dip or expose to liquid

Do not use the Monitor if there is any sign of damage

Use only PH Neutral Cleaning products.

This product is not suitable for mechanical re-processing or sterilization.

Monitor Cleaning Instructions: Only the Carry Case and if necessary the Monitor surfaces may be cleaned and/or disinfected. Use moist (not dripping) wipes with 70% solution of isopropyl alcohol, or very dilute Chlor-clean (1000ppm) or Chlorohexidine (1000ppm) or mild detergent, then allow to air dry naturally.

2.)Cleaning the SpO₂ probe.

Care:

Do not sterilize by high pressure, autoclave or washer

Do not dip the probe into liquid.

Do not use the probe if there is any sign of damage.

Use only PH Neutral Cleaning products.

This product is not suitable for mechanical re-processing or sterilization.

Cleaning instructions:

Use moist (not dripping) wipes with 70% solution of isopropyl alcohol, or very dilute

Chlor-clean (1000ppm) or Chlorohexidine (1000ppm) or mild detergent, then allow to air dry naturally.

3.) Cleaning of filter/water trap



Only the filter surface may be cleaned and/or disinfected. Use moist (not dripping) wipes with 70% solution of isopropyl alcohol, or very dilute Chlor-clean (1000ppm) or Chlorohexidine (1000ppm) or mild detergent, then allow to air dry naturally.

Replace the filter/water trap if it is dirty, shows any sign of contamination or if in any doubt about its condition.

7 Trouble Shooting Analysis

Simple analysis of problems

No.	Phenomena	Causes	Solution
1	The values of CO ₂ is reading too low, or 'OCCLUSION' appears on the screen.	<ol style="list-style-type: none"> 1. Leaking of filter or sampling tube 2. Occlusion of filter or sampling line 3. Out of Calibration 	<ol style="list-style-type: none"> 1. Check and replace filter or sample line 2. Clear the gas loop Occlusion 3. Re-calibrate using standard gas.
2	The values of CO ₂ is zero <ol style="list-style-type: none"> 1. Screen indicating PUMP ERR and big noise. 2. Screen indicating IR-LAMP-BAD 3. Screen indicating CO₂ SENSOR ERR 	<ol style="list-style-type: none"> 1. Internal leaking inside the Pump gas loop 2. The IR lamp resource of sensor damaged 3. IR Sensor broken 	Contact the Distributor or manufacturer for repair.
3	Screen indicating CAL-ERR	The last calibration has failed.	Re-calibrate using standard gas.
4	Screen indicating POWER-ERR	Damaged or incorrect power supply.	Contact Distributor or manufacturer.
5	The CO ₂ wave is not normal. <ol style="list-style-type: none"> 1. Screen indicating TEMP-HIGH 2. Screen indicating TEMP-LOW 3. Screen indicating TEMP-IMBALANCE 	<ol style="list-style-type: none"> 1. Temperature too high. 2. Temperature too low. 3. Sharp ambient Temperature change 	Use in normal environmental temperature range
6	No values of SpO ₂ or no wave	<ol style="list-style-type: none"> 1. Finger too cold 2. Interference of very strong external light 3. The measurement test of SpO₂ and blood 	<ol style="list-style-type: none"> 1. Warm up finger 2. Avoid strong external light. 3. Place SpO₂ sensor on other arm or

		<p>pressure are done on the same arm.</p> <p>4. Red light in the sensor no flashing.</p> <p>5. Infrared and collector of sensor is not clean</p>	<p>position.</p> <p>4. Renew SpO₂ sensor</p> <p>5. Clean internal parts of SpO₂ Sensor</p>
7	Flashing red color  and closed down automatically.	1. No Battery Charge.	1. Connect to Battery Charger.
8	Still flashing red color  after the power is supplied and AC indicator no light.	1. Battery Charger power working abnormally.	1. Check battery charger and cable and replace as necessary.

Attention: Please contact your distributor if you require advice, replacement parts and/ or service.

Appendix 1. Explanations of Terms in this Manual

MENU	Menu
EtCO ₂	The CO ₂ concentration of expiration end phase
INCO ₂	The CO ₂ concentration of inspiration phase
SPO ₂	Oxygen saturation
RR	Respiration rate
PR	Pulse rate
mmHg	Millimeters Mercury
kPa	Kilopascal
ALARM-H	Alarm high limit
ALARM-L	Alarm low limit
LINE	Line curve
FILL	Filled or solid under waveform
BEEP_VOLUME	Pulse volume
ALARM_VOLUME	Alarm volume
APNEA	Apnea or breathing stopped for a set period of time
BPM	Breaths per minute
SET	Setup
N ₂ O:	Nitrous oxide
HELIUM	Helium gas
O ₂ CONCENTRATION	O ₂ concentration compensation
ANAESTHETIC GAS	Anaesthetic gas
ZERO GAS	Base point or Zero point
BTPS	Temperature and deep lung air pressure compensation
CALIBRATE	Calibration
CANCEL:	Cancellation
OCCCLUSION	Blocked filter/water trap or gas sample line

Appendix 2. Changing compensation of balance gas

Attention:

Only the trained personnel may carry out the following the procedure. Contact your Supplier for training and advice.

Enter the engineer menu as follows:

Press **+** and **▼** two buttons simultaneously to enter the following menu.




Figure A2.1

In this menu, press **▲** or **▼** button to move the cursor up or down, press **+**

button or **-** button to change the data highlighted by the cursor.

Some items of this menu can be directly adjusted, such as LOAD-DEFAULT or EXIT: to press ENTER button, exit without saving or changing data. In this menu, press MENU button, then to exit this menu and enter the main menu.

This menu includes the following setups:

BALANCE GAS: AIR, N₂O, and HELIUM

O₂ CONCENTRATION: 20%-99%

ANAESTHETIC GAS: 0-20%

ZERO GAS: AIR, N₂

BTPS: ENABLE, DISABLE

MENU: UNLOCK, LOCK

LOAD DEFAULTS

CALIBRATE

Attention:

- 1.) When the menu is locked, this menu is disabled. To unlock the menu, press **+** and **▼** to enter engineer menu and change “unlock” to “lock” in the MENU setting. This is to avoid the misoperation of the patient against the preset of the doctor.
- 2.) CALIBRATE is for CO₂ concentration recalibration. Long press ENTER button for 8 seconds to enter this menu.

Default values are as follows:

BALANCE GAS: AIR

O₂ CONCENTRATION: 20 %

ANAESTHETIC GAS: 0 %

ZERO GAS: AIR

BTPS: DISABLE

MENU: UNLOCK

Appendix 3. Calibration of EtCO₂ Accuracy

Attention: Only trained personnel are allowed to carry out the following procedure. Contact your Supplier for training and advice.

The monitor has been calibrated before being shipped by the manufacturer. Generally the user does not need to calibrate this device other than the recommended annual check. To check the unit using Cal Gas the following procedure must be obeyed.

1. Required Parts and Items:

1.) **CO₂ standard gas** - Concentration is normally 5-8%

2.) **Three-way connector:** A three way connector with an inner diameter of 1-3 mm (one connection vented to open air) must be used to protect the monitor when calibrating using a CO₂ standard gas bottle see below figure . The device **will be damaged** by the high pressure of the standard Cal Gas Bottle if the connector is not used. It is strictly forbidden to connect the cal gas bottle directly to the device. One end of three-pass connector must be directly open to air to release gas pressure and protect the monitor.

3.) **Two tubes (whose length can extend outside room):** The standard gas flows into the air continuously through the three way connector and the module pump also vents the gas that is checked. During calibration CO₂ gas of a higher concentration can easily and quickly accumulate around the device. To prevent any potential of this affecting and influencing the calibration of the Zero base vent the connections from the three way adapter and the monitor to outdoor.

2. Connect as follows:




Figure A3.1

3. Warm-up

Turn on power and run the unit for 20-30 minutes and adjust the pump flow rate to over 120cc/min. To check if there is a leak use the following method: Squeeze the sampling tube by hand, the operating noise of the sampling pump will increase noticeably. If the sampling pump does not accelerate and its operating noise also does not change then there must be a leak in the gas loop. You must then find out where the leak is and solve it, otherwise, it will lead to incorrect calibration. After warming-up, open the flow of standard cal gas, and listen if the sound of pump is as same as original one. If the pump's turning is slow and its turning sound is weak, that means the standard gas pressure/flow is too large.

Turn down the Cal Gas flow rate until the sound of sampling pump resumes its original volume.

4. Calibrate

Enter the engineer menu (procedure given at Appendix 2,), highlight CALIBRATE, long press ENTER button for 8 seconds to enter the next menu.




Figure A3.2

Highlight STANDARD GAS and adjust the value to that of the concentration of CO₂ standard gas. If the standard gas concentration precision is to 2 decimal places numbers round up accordingly.

Then highlight CAL-BEGIN and long press the ENTER button for 8 seconds, at the same time, open the standard gas and the device will begin to calibrate. The screen will display the message 'ADJUSTING' as shown in below figure.




Figure A3.3

The thick cross bars in the display will be erased as time passes and the calibration will end when they are completely erased. If the calibration is successful, the menu will show ADJUST OK and subsequently exit into the main menu. If the calibration is unsuccessful, this menu will show ADJUST ERR. If this occurs the loop

needs to be checked to determine if there is a leak or standard gas has run out (the pressure indicator of gas bottle shows 0). The Calibration menu will remain if the calibration is unsuccessful.

If you require to exit this menu during the calibration press the MENU button or highlight CANCEL and press the ENTER button.

Note: Remember to close the valve of the standard gas to prevent wastage.

Appendix 4. Part Numbers and Consumables listing

CR-ASK900B Intubated Adult/Paediatric Airway Sampling Kit s x 10 Pack - each kit includes 1 of each item:

Water Trap/Filter T3 for PC-900B Monitor

Circuit Adaptor, 22F/15M with Gas Sample Port

Gas Sampling Tee, Male/Male Luers, 1.27mm ID x 3.0m

CR2500-0000218	Disposable Water Trap/Filter T4F (Female Luer Lock Connector)	Pack of 10
CR2500-0000240	Disposable Water Trap/Filter T4M (Male Luer Lock Connector)	Pack of 10
WL99370010	Elbow Connector Sampling Tee , 22F/15M with Port & and Cap, Adult/Paed	Pack of 50
QPS8003-50	Gas Sampling Line, Male/Male Luers, 1.27mm ID x 3.0m	Pack of 50
QO51035	Straight Connector Sampling Tee, 22F/15M with Port and Cap	Pack of 100
QO12090	Male to Male Luer Lock Connector to convert Water Trap/Filter T3 to Male Luer – allows use of Female Luer ended Sample Lines	Pack of 10

PLEASE NOTE: You will need to adjust the Flowrate of the sample pump down to 50 ml/min to allow the use of narrow bore very low flow Samples lines. Failure to do so may result in the Occlusion Alarm becoming active and/or premature Pump failure due to high resistance over stressing.

PB-331010	PRO-Breathe® CO2 Sampling Mask with O2 Delivery, Adult with 2.1m O2 Tubing, Female Luer connector	Pack of 50
MA4000	Nasal CO2 Sample Line, Adult, 2.1m with Male Luer	Pack of 25
MA4100	Nasal CO2 Sample Line, Paediatric, 2.1m with Male Luer	Pack of 25
MA4707	Nasal CO2 Sample Line with O2 Delivery, Adult, 2.1m with Male Luer	Pack of 25
MA4703	Nasal CO2 Sample Line with O2 Delivery, Paediatric, 2.1m with Male Luer	Pack of 25
CR15040050	SpO2 Sensor (Sub-D), Silicone, Adult, 2m	Pack of 1
CR15040051	SpO2 Sensor (Sub-D), Silicone, Paediatric, 2m	Pack of 1
CR15040022	SpO2 Sensor (Sub-D), Finger Clip, Adult, 2m Cable	Pack of 1
CR15040055	SpO2 Sensor (Sub-D), Finger Clip, Paediatric, 2m Cable	Pack of 1
CR15040017	SpO2 Sensor (Sub-D), Y Type with Silicone Wrap, 0.9m Cable	Pack of 1

CR2302-0000013	Replacement Lithium Battery for PC-900B Capnograph	Pack of 1
CR2903-2000010	Charger Cable (USB to mini USB) 1.5m	Pack of 1
ACA-USB2UK	5V DC Mini USB to UK Plug Adapter for use with above	Pack of 1
CR-12VDC9B	12V DC Vehicle Power Adapter to 5V DC Mini USB 3m length	Pack of 1
PROBAG-CAP	Heavy Duty Cushioned Carry Case for PC-900B Monitor	Pack of 1

Warning:

PLEASE USE ONLY GENUINE RECOMMENDED SPARE PARTS AND ACCESSORIES
OTHERWISE YOUR WARRANTY WILL BE INVALIDATED

Attention:

Please contact **your distributor if you require advice, replacement parts and/or service.**

Appendix 5. Guidance and manufacturer's declaration -

Electromagnetic compatibility

Table 1 Guidance and manufacturer's declaration-electromagnetic emission-for all
EQUIPMENT AND SYSTEMS

<p>This device is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment or system should assure that it is used in such an environment.</p>		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	<p>This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.</p>
RF emissions CISPR 11	Class A	<p>This device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.</p>
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC61000-3-3	Complies	

Table 2 Guidance and manufacturer's declaration-electromagnetic immunity for
all EQUIPMENT AND SYSTEMS

This device is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment or system should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
<i>Electrostatic discharge (ESD)</i> <i>IEC61000-4-2</i>	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	<i>Floors should be wood, concrete or ceramic tile. if floors are covered with synthetic material, the relative humidity should be at least 30%</i>
Electrical fast transient/burst IEC61000-4-4	± 2 kV for power Supply lines ± 1 kV for input/output lines	± 2 kV for power Supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line (s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 0,5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles	$< 5\% U_T$ ($> 95\%$ dip in U_T) for 0,5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the

<p>input lines IEC61000-4-11</p>	<p>70% U_T (30% dip in U_T) for 25 cycles</p> <p>5% U_T (>95% dip in U_T) for 5 s</p>	<p>70% U_T (30% dip in U_T) for 25 cycles</p> <p>5% U_T (>95% dip in U_T) for 5 s</p>	<p>equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.</p>
<p>Power frequency (50Hz/60Hz) magnetic field IEC61000-4-8</p>	<p>3A/m</p>	<p>3A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>
<p>NOTE: U_T is the a.c. mains voltage prior to application of the test level.</p>			

Table 3 Guidance and manufacturer’s declaration – electromagnetic immunity-for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an electromagnetic environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \text{ 80MHz to 800MHz}$ $d = 2.3\sqrt{P} \text{ 800MHz to 2.5 GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ,a should be less than the compliance level in each frequency range .b</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	















			<p>Interference may occur in the vicinity of equipment marked with the following symbol.</p> 
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p>			
<p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a: Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which this device is used exceeds the applicable RF compliance level above, this device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this device.</p> <p>b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.</p>			

Table 4 Recommended separation distances between portable and mobile RF communications equipment and the equipment or system-for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

<p>This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the equipment or system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the equipment or system as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150kHz to 80MHz $d = 1.2\sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2,5GHz $d = 2.3\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

	Manufacturer
	Date of manufacture
	Lot number
	Follow instructions for use
	Keep in a cool, dry place
	Imported by
	Product code
	Authorized representative in the European community
	WEEE disposal
	This item is compliant with Directive 93/42/EEC
	Caution: read instructions (warnings) carefully
	Keep away from sunlight

IPX1	Covering Protection rate
	Expiration date

GIMA WARRANTY TERMS:

The Gima 12-month standard B2B warranty applies.



Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.