

Scala

Polygraphy for varied uses.

Scala is primarily a standard **polygraphy device** for the detection of Sleep-Disordered Breathing that fulfills the requirements of the German EBM code for billable medical services. Furthermore, within the standard delivery Scala offers many possible uses at no extra cost for treatment of different user groups – all in one device. Given its variety of options, Scala is highly suited for use in a variety of fields, including pneumology, cardiology, ENT, sleep medicine and neurology.

Special features at a glance

Thermistor – In addition to a nasal cannula, a thermistor can be used to measure the patient's respiratory flow.

Snore microphone – A sensor to record snoring sounds can be applied to the larynx. It is possible to detect snoring at the same time by means of the dynamic pressure signal from the nasal cannula.

Differential pressure measurement – The integrated differential pressure measurement to manage PAP and NIV therapies offers very high-quality recording of the flow signal at high and quickly changing pressures.

Configurable electrode ports – The device has equipped with four electrode connection ports which can be freely configured and used for different electrophysiological tracings such as EMG, ECG or EEG.

Recording of respiratory movements – In addition to pressure pad technology, it is possible to record respiratory movements with induction belts (RIP) simultaneously or alternatively.

Removable Battery – The removable battery allows a quick hand-off of the device to the next patient without having to check the charging times.

Online and Offline – It is possible to carry out an online measurement instead of the standard offline measurement. An optional WLAN variant of Scala is available for data transmission from the patient recorder to a computer. The device can be operated with or without cables. A time-synchronized recording of video and audio signals can be conducted, as can a time-synchronized feed of external analog signals from therapy devices or transcutaneous CO₂ monitor.

No maintenance or service required – The device requires no maintenance or service. However, Löwenstein offers a regularly scheduled check of devices in order to ensure uninterrupted and trouble-free operations in daily use.



Scala



Device Description

Scala (standard delivery)

Article Nr.

1930041

Accessories

	Article Nr.		Article Nr.
Nasal Cannulae		Effort sensors	
Nasal cannula, 50 cm, luer lock	200-0312/10	Effort sensor short, adult	930393
Oro-nasal cannula, 60 cm, luer lock, small lip tab	500003	Effort sensor long, adult	930394
Oro-nasal cannula, 60 cm, luer lock, large lip tab	200-0312/15	Effort belt	
Sets		Size S (98 cm), blue	930260
Thermistor Set	1930320	Size M (122 cm), black	930261
Rechargeable battery Set	1930180	Size L (142 cm), red	930262
Snore microphone Set	1930321	USB connection cable, 1.8 m	930224
LEG Set	1930410	Pressure Diff Adapter Scala	930322
ECG Set	1930408	Set of double tubes for Pressure Diff Adapter, 50 cm	WM 96503
EEG Set	1930018	Flexible tube adapter for Pressure Diff Adapter	WM 95263
Patient bag	930253	SpO ₂ Softclip finger sensor	930300
Scala mount	930250	Hook-and-loop armband	200-0311
Rechargeable battery for Scala	930180	Transport case	927304
Rechargeable battery cover for Scala	930142	Scala cleaning set	930376
Charging tray for Scala	930204		
Power supply unit for battery charger for Scala	930209		

Technical data

Dimensions W x H x D:	135 mm x 30 mm x 77 mm
Weight:	230 g, including rechargeable battery
Storage medium:	Internal storage
Storage capacity:	Min. 100 hours
Error indicator:	LED on top of device
Power supply:	Rechargeable Li-Ion battery with integrated charging electronics
Data transmission:	USB interface with cable for data transmission
Sampling rate:	1 Hz to 1000 Hz with 24-bit resolution
Length of recording:	Up to 24 hours (depending on recording configuration)
Maintenance and service:	No maintenance or service required
Software updates:	Free software updates for our analysis software for the life of device
Software license:	No software license required

CE 0483



Dr. Fenyves und Gut Deutschland GmbH
Daimlerstr. 23
72414 Rangendingen, Germany

Sales + Service
Löwenstein Medical
Arzbacher Straße 80
56130 Bad Ems, Germany
T. +49 2603 9600-0
F. +49 2603 9600-50
info@loewensteinmedical.com
loewensteinmedical.com

