

bio roller® G10

ALL IN ONE KIT
Microneedle rolling system



Advantages

| | | | | |
|--|---|--|--|--|
| <p>720 Needles 1.0 mm</p> <p>For skin care on face</p> | <p>300 Needles 0.5 mm</p> <p>For skin care on eye contour</p> | <p>27 Needles 1.5 mm</p> <p>For skin care on topical</p> | <p>Ice roller head</p> <p>For skin ice massage</p> | <p>Hedgehog roller ball</p> <p>For skin massage</p> |
| <p>Double jade roller</p> <p>For skin massage</p> | <p>Steel ball</p> <p>For skin massage on topical</p> | <p>Silicone brush</p> <p>For skin cleaning</p> | <p>Essence spoon</p> <p>For applying essence to skin</p> | <p>Disinfection basin</p> <p>For sterilizing the roller head</p> |

G10 is equipped with 10 different accessories that can be matched in different combinations for different skin problems. The unique curved surface massage roller system can massage the skin well and deeply soothe the skin. Different accessories of the roller microneedle can bring different effects, and can be used on the scalp, face, and various parts of the body, including sensitive skin around the eyes, local area, and neck. It can solve your skin problems.


DIFFERENT ROLLERS ARE USED FOR DIFFERENT FUNCTIONS

- Shrink pores, Wake up the skin** (Ice roller head)
- Cheek massage, Jade skin care** (Double jade roller)
- Massage acupoints, Metal press** (Steel ball)
- Eliminate neck line, Massage body** (Hedgehog roller ball)
- Nose, Lip groove, Mid-brow line** (720 Needles 1.0 mm)
- Eye brow line, Crow's feet** (300 Needles 0.5 mm)
- Cheek, Lips** (27 Needles 1.5 mm)

Easy to assemble



Non-slip handle



Disinfection tank




Better Skin Care




EFFECTS OF USING MASSAGE ROLLERS

There are 4 different massage roller heads in G10 kit each with a different effect. You can care and massage every inch of skin to achieve the effects of awakening the skin, eliminating puffiness, shrinking pores, relaxing muscles, eliminating puffiness, shrinking pores, and calming the skin. For puffy faces, you can also apply the serum to the skin with a massage roller.



| | | |
|-------------------|-------------------|------------------|
| | | |
| Wake up your skin | Stamp and massage | Massage the body |




EFFICIENT SYSTEM FOR ANTIAGING SKIN REGENERATION

A derma roller contains hundreds of microneedles. Micro-injuries are created by the microneedles thus causing collagen and elastin to grow, resulting in skin tightening. Microneedling is also effective in treating acne scars and hyper-pigmentation.



| | |
|---------------------|--------------------|
| | |
| Reduce eye wrinkles | Reduce crow's feet |

Product Parameters



Specification

| | | | |
|--|---|------------------------------|---|
| Brand | Bio Roller | Model | G10 |
| G.W. | 17.3*10.3*4.2 cm/1 pc/ 340 g | 1 carton | 54 * 37.5 * 46 cm/100 pcs/36.5 kg |
| Double jade roller material : Natural jade stone | | Needles material : Stainless | Face Brush material : Silicone |
| Ice roller head material : Metal | | Handle Material : ABS | Disinfection basin Material : ABS |
| Hedgehog roller ball Material : ABS | | Steel ball material : Metal | Essence spoon material : ABS |
| Customization (OEM & ODM) | logo private label package new shape change color | Packing List | 1. Box 2. Tray 3. Manual 4. Handle 5. 720 Needles 1.0 mm 6. Ice roller head 7. Hedgehog roller ball 8. Double jade roller 9. Silicone brush 10. 300 Needles 0.5 mm 11. 27 Needles 1.5 mm 12. Steel ball 13. Essence spoon 14. Disinfection basin |