

120 Needles

Adjustable Needle Size

Medical grade stainless steel needle
Needles arranged neatly
Needle Straight & Sharp



Bio needle 120 punctures the skin and creates thousands of microchannels in a short period of time. Used in conjunction with skin care products and hair growth essence products, the essence and active ingredients are delivered to the deep layers of the skin to enhance absorption, whiten the skin, awaken hair follicles, and promote hair growth.

Advantages



Adjustable needle length for easy operation

Can be applied to different skin areas

0-2.0 mm Rotate to adjust the needle length

0 0.5 1.0 1.5 2.0

Turn counterclockwise to change the length from short to long




Sharp needle tip Evenly distributed needles

Quickly open the skin's absorption channels
Allows active ingredients to effectively penetrate into the skin



Medical grade steel needle, light texture and high strength

Blister Packaging



Keep needles safe during transportation.

Protective cap



Prevent bacterial invasion, safe and hygienic.

Grow hair

Open the micro-absorption channels of the scalp

Help delivery and absorption of hair growth serum, improve traditional hair growth.


Arc edge design, more suitable for the skin



Bio needle 120 stimulates the skin and creates thousands of microchannels in a short period of time. Used in conjunction with hair growth essence, it transports growth factors to the deep layers of the scalp, enhances blood circulation in the scalp, activates damaged hair follicles, and promotes hair growth from the roots.


It's not that hair growth serum has no effect But the introduce method is not ideal

bio needle 120 penetrates the surface of the scalp, allowing the nutritional essence to truly penetrate



The absorption rate of regular application is 15%

Regular application of essence acts on the surface of the scalp and is difficult to be absorbed.







The absorption rate of microneedle introduction is 95%

Micro-needle opens the absorption channel, allowing the essence to penetrate into the deep layer of the scalp.

Before and After

Before

After





Note: bio needle 120 is only an auxiliary tool for essence introduction. The actual effect varies from person to person, and the final effect will also vary depending on the hair growth essence used.

Skin Care


Solve various skin problems

Break through the stratum corneum and introduce skin care essence, allowing the active ingredients to deliver deep into the skin.


- Reduce stretch marks
- Repair sensitive skin
- Reduce deep scars
- Reduce acne scars
- Introduce essence
- Reduce fine lines
- Improve rosacea
- Reduce pigment
- Reduce wrinkles
- Firm the skin
- Smooth skin
- Lift the skin



Red bloodshot



Reduce acne marks



Reduce scars

Different lengths for different skin problems

microneedle penetrates deeply to improve skin care effects



0.25-0.5mm

Red bloodshot, sensitive skin, fine lines around eye contour

0.5-1.0mm

Acne marks, fine lines, large pores

1.0-1.5mm

Acne pits, pigment, scars

1.5-2.0mm

Reduce stretch marks

2.0mm

deep scars

*The thickness of facial skin is about 1.5 mm. Needles longer than 1.5 mm are not suitable for use on facial skin.

Product Parameters

Specification

Brand	Bio Needle	Model	120 Needle
Needle Adjustment Range	0-2.0mm	Needle Number	120 Pins
Handle Material	ABS	G.W.	150*53*30 cm / 1 pc / ? g
Needle Material	Stainless steel	1 Carton	31.5*28*31.8 cm / 100 pcs / ? kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box ② 120 Needle ③ Inductive base

