

**biopen i q2**  
 Skin Care Instrument



IQ2 is a beauty instrument that integrates microneedling, EMS electrostimulation, and LED light therapy in a single device for beauty salon and home use. It is based on an advanced laser, it uses smart chips to adjust energy output to the skin, to protect the user's skin and make the treatment more effective and safer.

**Focus on microneedling therapy in beauty and cosmetics industry**  
 Facial introduction, auxiliary absorption

- Over 10 years Professional R&D team
- 100,000+ cases and free workshop
- Brand name and original manufacturer of well-known brand (Derm)®



**QUALIFICATION CERTIFICATE**

- ISO9001:2015 Quality Management System
- CE (EMC) Certification
- CE (LVD) Certification
- CE (RoHS) Certification
- CE (REACH) Certification
- CE (FDA) Certification
- CE (FDA) Certification
- CE (FDA) Certification
- CE (FDA) Certification

**Selling point**


**Triple effects Reluvenate the skin**

Microneedling EMS Electroperoration LED skin care



**25 Pole Pulse EMS Electroperoration**  
 Stable and average energy  
 Three-dimensional stimulation


**EMS electroperoration**  
 EMS needle  
 Deeply penetrates and connects  
 microcurrent energy




**EMS Microcurrent energy field**  
 EMS microcurrent energy field  
 EMS microcurrent energy field  
 EMS microcurrent energy field

**Four skincare in one use**  
 Microneedling-EMS Electroperoration-Blue Light-Red Light

Support skin by  
 Microneedling EMS Electroperoration 430nm Blue Light 630nm Red Light



**For Skin Care**



- ✓ Fine lines, wrinkles
- ✓ Pore size, acne scars
- ✓ Skin rejuvenation
- ✓ Improve skin texture
- ✓ Freshen up on daily
- ✓ Large vessel closure
- ✓ Smooth neck on thigh

**For Hair Growth & Beard Growth**

It's not that the hair restorer has no effect, but your deliver method is not ideal. Say goodbye to traditional skin care methods.

Good the hair restorer and also not good the deliver method is not ideal. Say goodbye to traditional skin care methods.



**Microcurrent absorption**  
 Microcurrent energy field  
 Microcurrent energy field

**Advantages**

Power when it is placed and starts when off

Smart touch screen



LED has better working to help skin from the energy battery power. Microcurrent energy field is better than the traditional.

Low battery warning



**Product Parameters**



Model	Power	Material	Color	Weight
Biopen i q2	10W	Aluminum	White	150g
Biopen i q1	10W	Aluminum	White	150g
Biopen i q3	10W	Aluminum	White	150g
Biopen i q4	10W	Aluminum	White	150g
Biopen i q5	10W	Aluminum	White	150g
Biopen i q6	10W	Aluminum	White	150g
Biopen i q7	10W	Aluminum	White	150g
Biopen i q8	10W	Aluminum	White	150g
Biopen i q9	10W	Aluminum	White	150g
Biopen i q10	10W	Aluminum	White	150g