


SE-2003/SE-2012

Holter System



EDAN

A world of potential

SE-2003/SE-2012 Holter System

EDAN holter system offers you a lightweight and compact digital recorder, as well as fast and powerful analysis software. Optimized workflow and adjustable parameters aim to meet your different clinical needs.

Holter Recorder

- Lightweight design: only 42g
- 3/12-channel digital recording
- Real-time waveforms display
- Data stored in SD card
- Pacemaker detection
- Up to 7 days recording with only one AAA battery



Holter Software

Data Processing

- Full disclosure of 24/48/96 hours and 7 days ECG data
- Download data via SD card or USB cable
- Automatic and manual edit of patient data
- Graphics and trends display

EMR Compatibility

- GDT/XML protocol support

Powerful Analysis

- QRS template classification
- Arrhythmia analysis
- ST segment analysis and myocardial ischemia evaluation
- Atrial fibrillation analysis
- Obstructive sleep apnea syndrome analysis
- Heart rate turbulence (HRT) analysis

