

# CardioChek® Plus

## Fast, Accurate Lipid and Glucose Results at the Point of Care


Recognizing the strong connection between heart attack, stroke, and diabetes, the CardioChek® Plus analyzer provides rapid on-site blood chemistry values and measurements to help identify at-risk individuals. It's one of the fastest, most cost-efficient, and user-friendly analyzers. The CardioChek Plus system accurately tests lipid profile and glucose simultaneously with one fingerstick. In as little as 90 seconds, it provides on-site results using the same technology as clinical laboratories.

### Measures:

- Total Cholesterol
- HDL Cholesterol
- Triglycerides
- Glucose

### Calculates:

- LDL Cholesterol
- TC/HDL Ratio
- LDL/HDL Ratio
- Non-HDL Cholesterol



Simultaneous testing of lipid profile and glucose

### Fast

On-site lipid and glucose results in as little as 90 seconds

### Portable

Compact and lightweight, easy to transport and store  
Handheld, battery-operated

### Patient-Friendly

Fingerstick with small sample

### Wireless Communication

Easy data transfer

### High ROI

Rapid throughput and low waste, with minimal start-up cost  
Fully reimbursable



### Certified Accuracy

CLIA-waived

Meets NCEP guidelines for accuracy and precision

FCC and ISO certified wireless communications



# Wide Variety of Test Strips Provides Medical Professionals Flexibility

## The CardioChek Plus Lipid+eGLU™ SmartBundle™ includes:

- (3) 15-count Lipid Panel and
- (1) 50-count eGLU (electrochemical glucose) test strip vials.



Three-Analyte



### Lipid Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Triglycerides  
Calculates: LDL Cholesterol, Total Cholesterol/HDL Ratio, Non-HDL Cholesterol



### CHOL+HDL+GLU Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Glucose  
Calculates: Total Cholesterol/HDL Cholesterol Ratio



### METABOLIC CHEMISTRY Panel Test Strip

Measures: HDL Cholesterol, Triglycerides & Glucose

Dual-Analyte



### CHOL+GLU Panel Test Strip

Measures: Total Cholesterol & Glucose



### CHOL+HDL Panel Test Strip

Measures: Total Cholesterol & HDL Cholesterol  
Calculates: Total Cholesterol/HDL Cholesterol Ratio

Single-Analyte



### TOTAL CHOLESTEROL Test Strip

Measures: Total Cholesterol



### HDL CHOLESTEROL Test Strip

Measures: HDL Cholesterol



### TRIGLYCERIDES Test Strip

Measures: Triglycerides



### GLUCOSE Test Strip

Measures: Glucose



### eGLU Test Strip

Measures: Electrochemical Glucose

Single-analyte test strips available as special order include Creatinine and Ketone.



**PTS Diagnostics**  
7736 Zionsville Road, Indianapolis, IN 46268 USA  
Office: +1 (317) 870.5610 FAX: +1 (317) 870.5608  
Toll-Free in the U.S. : +1 (877) 870.5608

[ptsdiagnostics.com](http://ptsdiagnostics.com)