



SE 315  
SE 75      SE 325  
SE 85      SE 335  
SE 95



SE 345      SE 355      SE 365



45620

### Illuminated Test Types

45620 Internally full illuminated test cabinet with 2 concealed tubular light fittings each 230V 60W (no. 45650) and cable, with two test type panels by choice.

- SE 315: numerals for 5 m
  - SE 75: mirror writing
  - SE 325: numerals for 5 m
  - SE 85: mirror writing
  - SE 335: letters for 5 m
  - SE 95: mirror writing
  - SE 345: Landolt ring for 5 m
  - SE 355: E-hooks for 5 m
  - SE 365: pictures for children 5 m
- Dim.: 410 x 560 x 105 mm (W x D x H);  
weight: 4.2 kg

45650 Spare light tube 230V 60W for 45600/45620

45401 Optotype chart, single

42700 Red/green glasses (not shown in the catalogue)

42710 Red/green in hand-held frame (not shown in the catalogue)

### Colour Vision Tests

- For colour defects
- In four different types
- Book size

47640 Colour vision tests by Ishihara, 38 plates

47600 Colour vision tests by Ishihara, 24 plates

47610 Colour vision tests by Ishihara, 14 plates

47620 Colour vision tests by Ishihara, for illiterates, 10 plates


- For examination of colour defects
- 10 charts in book size
- Already available for preschool age

47625 Colour vision tests by Matsubara for children


Other colour vision tests

47650 Colour test panel 15 D by Farnsworth (saturated)

47651 Farnsworth colour test 15 HUE (desaturated)



47600



47625

### Oculus Illuminated Test Types

45600 Illumination cabinet with two concealed tubular light fittings no. 45650 ea. 60W, 1 test chart as per choice (46016 to 46926). Only 230V.

Size: 410 x 560 x 105 mm (W x D x H); weight: 2.5kg



45600 with  
46346

### Visual Acuity Charts For Distance

The version of the visual acuity charts is designated by the last digit of the order no. (3 or 6). When ordering, please add the required last digit to complete the order no. V/A = visual acuity for 5 m test distance (the acuities for 6m test distance are also printed on the charts, except 4693: only 5m).

Version 3 = paper

Version 6 = plastic covered

Size approximately 290 x 490mm

4601. Letters, available in versions 3, 6; V/A 0.125-1.0

4603. Pictures for children available in versions 3, 6; V 0.125-1.0

4606. Numbers, available in versions 3, 6; V 0.125-1.0

4611. Letters, available in versions 3, 6; V 0.1-1.25

4616. Numbers, available in versions 3, 6; V 0.1-1.25


4634. Numbers/letters, available in versions 3, 6; V 0.125-1.66

4692. Illiterate E's, available in versions 3, 6; V 0.1-1.0

4693. Pict. f. children by Loehlein, available in versions 3,6; V 0.2-1.0

When ordering, please specify the order number complete with the 5<sup>th</sup> digit. (Example: 4601 3 = Letters, paper, norm. writing).


Unless otherwise stated we supply version 6: normal print, plastic covered.




4601




4603




4606



4634



4692



4693



47300

### H-Test

- Eye test for children with musical feedback
- From the age of three, also for children with speech problems

47300 H-Test, child's visual acuity test, acc. to Hohmann/ Haase for children from the age of three with musical feedback

47301 Test forms for H-Test (5 x 100 sheets)

47304 Visual acuity chart for H-Test

- Qualified for examination for children until the age of six

### Bino-Near Vision Tester

**For the examination of all important near-vision functions of the two eyes together and of each eye separately**

- Measurement of the monocular and binocular visual acuity at a distance of 40 cm
- Measurement of the near point position
- Test for near fusion
- Test for near convergence
- Test for near stereoscopic vision
- Test for near accommodation balance

45900 Near Vision Tester incl. base for the Bino-Near Vision Tester and plug-in power supply unit (without test chart)

45911 Test chart with 6 tests

35080 Coloured lens, red

45902 Spare bulb 6V/3W



45900



45911

### Visual Acuity Test With Landolt Rings

**For near and distance  
(C-Test by Haase/Hohmann)**

Near vision test consisting of 6 charts with logarithmically stepped visual acuities (40cm reading distance). The feature of this vision test is: the spacing (2 different values) between the optotypes is maintained constant for all visual acuity levels, the smaller spacing value is comparable to text spacing. The direction of the C opening varies on a random basis in each line. The C diminishes in size from line to line and thus the visual effort required increases. This near vision test can also be performed with children who are not able to read.

It is practical to use the demonstrator with Landolt ring: Oculus Order No. 46970.

47020 Visual acuity test (C-test) for near with Landolt rings, set of 6 charts in ring binder, 50 recording charts

Distance sight test consisting of 3 separate charts (5m reading distance). The charts for distance visual acuity testing feature the same parameters as for near vision, acuities 0.16-1.25

46966 Visual acuity test for distance (C-test) set of 3 charts, washable, plastic covered.

Size: 320 x 640 mm

46967 C-test single, chart No. 1 - rows 2.  
6' optotype-spacing

46968 C-test single, chart No. 2 - rows 17.  
2' optotype-spacing

46969 C-test single, chart No. 3 - single Landolt rings

### Near Vision Test Types

47130 Near vision test, Nieden 1-13, Illiterate E's, Landolt rings, newsprint, timetable, notes.  
Matt finish in plastic folder.

47170 Near vision chart, single E's and Landolt rings  
laminated, single sheet

47211 Near vision test by Aulhorn  
laminated



47020



47130