



## > THERMOMETRY

### Electronic Thermometry

- > SureTemp Plus 692 \_\_\_\_\_ 96
- > SureTemp Plus 690 \_\_\_\_\_ 96

### Ear Thermometry

- > Braun ThermoScan PRO 4000 \_\_\_\_\_ 98

## 96 Electronic Thermometry

### THERMOMETRY



**01692-200**  
SureTemp® Plus 692  
Electronic Thermometer

### SureTemp® Plus 692

Combines the speed of ear thermometry with the accuracy and reliability of traditional oral thermometry. Takes approximately 4- to 6-second oral, 10-second pediatric axillary, 15-second adult axillary and 10-second rectal temperatures.

- > Ergonomic design—comfortable to hold, easy to use, and very durable
- > Interchangeable, removable probe wells reduce the risk of cross-contamination
- > Large backlit LCD display with location identification capability
- > Includes a 60-second pulse timer for counting pulses and respirations
- > Your choice of security alarm systems to help prevent theft
- > Takes approximately 5,000 readings on 3 AA batteries

**01692-200** SureTemp® Plus 692 Electronic Thermometer with Wall Mount, Security System with ID Location Field, and Oral Probe with Oral Probe Well

**01692-201** SureTemp® Plus 692 Electronic Thermometer with Wall Mount, Security System with ID Location Field, and Rectal Probe with Rectal Probe Well

**01692-300** SureTemp® Plus 692 Electronic Thermometer, One Per Room, with Wall Mount, 9 ft Oral Probe with Oral Probe Well

**01692-301** SureTemp® Plus 692 Electronic Thermometer, One Per Room, with Wall Mount, 9 ft Rectal Probe with Rectal Probe Well

**01692-700** SureTemp® Plus 692 Electronic Thermometer, with Wall Mount, 9 ft Oral Probe with Oral Probe Well, with Rolling Stand and Holder for Spare Probe and Well

### SureTemp® Plus 690

Combines the speed of ear thermometry with the accuracy and reliability of traditional oral thermometry. Takes approximately 4- to 6-second oral temperatures, 10-second pediatric axillary and rectal temperatures, and 15-second adult axillary temperatures.

- > Interchangeable, removable probe wells reduce the risk of cross-contamination
- > Large LCD displays temperatures in Fahrenheit or Celsius
- > Last temperature recall
- > Convenient storage for 25 probe covers
- > Takes approximately 6,000 readings on 3 AA batteries



**01690-400**  
SureTemp® Plus 690  
Electronic Thermometer,  
Wall Mount, Oral Probe  
with Oral Probe Well

**01690-200** SureTemp® Plus 690 Electronic Thermometer, Oral Probe with Oral Probe Well

**01690-201** SureTemp® Plus 690 Electronic Thermometer, Rectal Probe with Rectal Probe Well

**01690-300** SureTemp® Plus 690 Electronic Thermometer, One Per Room, Wall Mount, 9 ft Oral Probe with Oral Probe Well

**01690-301** SureTemp® Plus 690 Electronic Thermometer, One Per Room, Wall Mount, 9 ft Rectal Probe with Rectal Probe Well

**01690-400** SureTemp® Plus 690 Electronic Thermometer, Wall Mount, Oral Probe with Oral Probe Well

**01690-700** SureTemp® Plus 690 Electronic Thermometer, Wall Mount, 9 ft Oral Probe with Oral Probe Well, with Rolling Stand and Holder for Spare Probe and Well

## SureTemp® Plus 692 & 690 Accessories



**05031**  
Disposable  
Probe Covers

**05031-101** Disposable Probe Covers (1,000 covers, packaged 25 per box)

**05031-110** Disposable Probe Covers (10,000 covers, packaged 25 per box)

**21326-0000** Wall Holder for 690/692

**21330-0000** Holder for Spare Probe & Probe Well for 690/692

**06138-000** Calibration Key, Assembly for 690/692

**02892-000** Probe and Well Kit, 4 ft Rectal

**02892-100** Probe and Well Kit, 9 ft Rectal

**02893-000** Probe and Well Kit, 4 ft Oral

**02893-100** Probe and Well Kit, 9 ft Oral

**02891-0000** Probe Well, Blue for 690/692

**02891-1000** Probe Well, Red for 690/692

**08273-000** Thermometer Stand, 690/692 with Holder for Spare Probe and Well

**21333-0000** Bracket, M600 Stand Upgrade, 690/692 with Holder for Spare Probe and Well

**01802-110** 9600 Plus Calibration Tester



**21330-0000**  
Holder for Spare  
Probe and Well



**02893-000**  
Oral Probe and  
Well Kit, 4 ft



**02891-0000**  
Oral Probe  
Well



**02891-1000**  
Rectal Probe  
Well

## 767 Wall System Add-on

**01767-200** SureTemp® 767 Wall Thermometer

**01767-600** SureTemp® 767 Wall Thermometer with Backlight



**01767-200**  
SureTemp® 767  
Wall Thermometer  
(shown with  
767 Integrated  
Diagnostic  
System)

# 98 Ear Thermometry

THERMOMETRY



**04000-200**  
Braun ThermoScan®  
PRO 4000 Ear Thermometer

## BRAUN *ThermoScan* PRO 4000

Provides the speed and accuracy you need for all types of patients. The exclusive ExacTemp™ technology provides an active user feedback system to ensure proper positioning and improved accuracy.

- > Improved accuracy and repeatability through a preheated probe tip and sensor
- > Probe cover detection system ensures precise infrared heat transfer for accurate temperatures
- > Automatic probe cover eject button for quick and easy disposal of used probe covers
- > Easy-to-read, backlit LCD displays temperatures in Fahrenheit or Celsius
- > Three-year limited warranty

**04000-200** Braun ThermoScan® PRO 4000 Ear Thermometer (without rechargeable batteries)

**04000-300** Braun ThermoScan® PRO 4000 Ear Thermometer with Recharging Base Station

**04000-400** Braun ThermoScan® PRO 4000 Ear Thermometer with Roll Stand

### Braun ThermoScan® PRO 4000 Accessories

**05075-800** Braun ThermoScan PRO 4000 Disposable Probe Covers (800 covers, packaged 20 per box)

**05075-005** Braun ThermoScan PRO 4000 Disposable Probe Covers (5,000 covers, packaged 20 per box)

**24001-0000** Braun ThermoScan PRO 4000 Recharging Base Station with Rolling Stand

**24001-1000** Braun ThermoScan PRO 4000 Recharging Base Station and Probe Cover Refill Dispenser

**53009-000** Braun ThermoScan PRO 4000 Alkaline Battery

**53020-0000** Braun ThermoScan PRO 4000 Rechargeable Battery Pack

**01802-110** 9600 Plus Calibration Tester

