

NEW ASKIR 30 is a desk-type electric suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. Suitable for tracheotomized patients, minor surgical applications and post-operative therapy in home-care.

# NEW ASKIR 30

Antibacterial and Hydrophobic Filter

Silicone Autoclavable Tubes

Vacuum Regulator

Vacuum Indicator (kPa and Bar)

Collection Jar in Autoclavable PC with Overflow Valve System



## TECHNICAL FEATURES

CE 0123

Optional Stand (ref. 27730)



Motor	Oiless and maintenance-free piston pump
Power Feeding	230V-50Hz Other configurations available upon request
ISO 10079-1 Classification	HIGH VACUUM / HIGH FLOW
Max Vacuum (adjustable)	-0.80 Bar -80 kPa -600 mmHg
Max free air flow rate	40 l/min
Noise Level	60,5 dB
Power consumption	107 VA
Fuse	1 x F1.6 A 250V
Duty cycle	Non-stop operation
Weight	Kg 3.60
Size	cm 35 x 21 x 18
Years of Warranty	2
Shipping carton	4
Place of Manufacturing	Italy

## MAIN APPLICATIONS

General Surgery

Endoscopy

Obstetrics

Veterinary

CA-MI

CA-MI srl Via Ugo La Malfa 31 - 43010 Pilastrò (PR) - ITALY  
Tel +39 0521 637133 - 631138 Fax +39 0521 639041  
export@ca-mi.it - vendite@ca-mi.it - [www.ca-mi.it](http://www.ca-mi.it)

SURGICAL  
SUCTION  
UNITS