

## ∴ HEINE BETA® 4 NT Rechargeable Handles and HEINE® NT 4 Table Charger

Next-generation energy management



∴ **BETA 4 NT** and **BETA 4 SLIM NT**  
for NT 4 table charger.

∴ **Li-ion technology:** no „memory effect“.  
Charge your handle how you like without  
worrying about charge level or battery condition.



∴ **High capacity:** more than double the capacity  
compared to previous NiMH rechargeable batteries.

∴ **High quality metal construction** for strength  
and durability.

∴ **Quick charge:** only approx. 2 hours to recharge  
with NT 4 table charger. Charge Status Indicator:  
light pulses from the illuminated charging wells  
indicating the charging mode. Continuous light  
indicates that the charging process has been  
completed.



∴ **Charge Status Indicator on  
the BETA 4 NT bottom insert:**  
the indicator illuminates in orange  
when the battery needs to be  
recharged.

BETA 4 NT rechargeable handle for NT 4 table charger X-007.99.396

NT 4 table charger X-002.99.494

BETA 4 SLIM NT rechargeable handle for NT 4 table charger X-007.99.474



## HEINE BETA® 4 USB Rechargeable Handle

Next-generation energy management



USB

• **Maximum safety and flexibility.** You can charge your HEINE BETA4 USB Handle with any USB plug or other USB power source, with absolute peace of mind. The integrated, exclusive Galvanic Separation Module protects the handle, and the user, even if a faulty or inferior plug is used, preventing the possibility of the full mains voltage flowing directly to the handle.

• **Li-ion technology:** no “memory effect”. Charge your handle how you like without worrying about charge level or battery condition.

**LI ION**

• **High capacity:** more than double the capacity compared to previous NiMH rechargeable batteries.

• **Quick charge:** approx. 4 hours to recharge.

• **High quality metal construction** for strength and durability.



• **Charge Status Indicator:** the indicator illuminates in orange when the battery needs to be recharged, blinks green during the charging process and turns permanently green as soon as it is fully charged.

• **BETA4 USB Li-ion rechargeable handle, incl. USB cord with medical approved plug-in power supply, handle rest** X-007.99.388

• **BETA4 USB Li-ion rechargeable handle (without USB cord)** X-007.99.387

