

SPHYGMOMANOMETERS WITH GAUGE ON THE CUFF

CE 0051

32724

Quality
linen
cuffTOP
QUALITY

Extra large dial Ø 65 mm

Extra large bulb
9x5 cm

• 32724 BERLIN SPHYGMO - velcro cuff

High quality sphygmomanometer with gauge mounted on cuff:

- high reliable manometer with extra large dial Ø 65 mm
 - top quality and resistant linen cuff with simulated leather support and velcro closure
 - extra large PVC bulb, 9x5 cm, latex free
- Made in Italy - German manometer

CE 0051

32723

ITALIAN
QUALITY

• 32723 NEW EXACTOPHON SPHYGMO - pull-on strap cuff

GIMA quality and reliability in a stylish and long-life sphygmomanometer. EXACTOPHON has a reliable manometer fixed on a simulated leather support, on a highly resistant linen pull-on straps cuff. Supplied in a nice carrying case and box. Made in Italy: GIMA.

- 32805 EXACTOPHON BLADDER - spare
- 32806 EXACTOPHON CUFF - spare

CE 0476



ECONOMIC

LATEX FREE

32719

• 32719 GIMATONO - latex free

Economic sphygmo with chrome-plated manometer fixed on the cuff, provided with a rubber protective ring. Nylon cuff with simulated-leather support. Two tubes adult size latex free TPU bladder, with standard latex free PVC bulb. Supplied in nylon bag and carton box. Made in P.R.C.