

Creative PC-300 Dual Mode Spot Check Monitor



Free Android Phone App

- Bluetooth enabled Multi Parameter Spot Checker
- Measures SpO2, NIBP, Pulse Rate and Temperature
- 4.3 inch LCD display
- Integrated Lithium battery
- Stores up to 6000 NIBP readings
- Able to measure temperature in 5 seconds
- Includes free Windows PC data management download software for Bluetooth and USB download
- Data upload via bluetooth for analysis and printout
- Optional ECG and Glucometer (available soon)
- Bluetooth Version 2.0
- 1 Year Warranty



Phone not included




PROTOCOLS AVAILABLE


SIGMA MEDICAL

Dedicated to Health

Part No.	NPC Code	Description	Qty	Price
CR-PC301	---	PC-301 NIBP Monitor (USB/Bluetooth v2.0 Data Transfer) with Adult BP Cuff, USB Download & Charging Cables	Single	£290.00
CR-PC302	---	PC-302 SpO2 & NIBP Monitor (USB/Bluetooth Data Transfer v2.0) with Adult Sensor, Adult BP Cuff, USB Download & Charging Cables	Single	£305.00
CR-PC303	---	PC-303 SpO2, NIBP & Temperature Monitor (USB/Bluetooth Data Transfer v2.0) with Adult Sensor, Adult BP Cuff, Tympanic Probe, USB Download & Charging Cables	Single	£325.00



Multiple parameters stored and recorded. ECG, SpO2, NIBP and Temperature.



Wireless information recorded in real time.



Plethymographic waveform information available for easy print outs.