

# LED MAGNIFYING LAMP - 5 DIOPTRES

## • 30737 LOUPE LED MAGNIFYING LIGHT

Smartly designed cold light LED lamp with 5-dioptres magnification. 5 wheels base and articulated arm allow easy movement. Manual in GB, FR, IT, ES, PT, DE  
Made in Europe.

**MIMSAL**



### TECHNICAL SPECIFICATIONS

Type of light:	LED
Dioptres:	5
Lens:	Ø 12.7 cm
Operating voltage:	230 V - 50 Hz
Class II	
Electrical security:	UNI EN 60601-1 (1990)
Electromagnetic compatibility:	EN 60601-1-2 (1993)



CE  
93/42

## • 30736 DERMALED MAGNIFYING LAMP

Cold light LED lamp with 5-dioptres magnifying lens and brightness adjustment. It comes with a high stability 5 wheels base and an articulated arm for easy movement and fixing of the lamp.  
Suggested for dermatology and beauty treatments.  
Made in China.



Height  
50+50 (arm) +  
95 cm (trolley)

**30737**



### TECHNICAL SPECIFICATIONS

Type of light:	LED
Voltage:	8 W
Dioptres:	5
Operating voltage:	100-240 V - 50/60 Hz
Net/Gross weight:	10-12 kg



**30736**