

## HEINE® mPack mini

Accurate examinations made easy – all day long.

The new, high performance power source **HEINE mPack mini** for LoupeLight 2 and MicroLight 2 offers freedom from cables that tie the users to desk, wall or static power sources. With a weight of only 81 g, this mobile power source is now approx. 70 % smaller and lighter than its predecessor.

The pouch which is included in the delivery can be attached to the belt and the mPack mini can also be charged while clipped to the belt. Thanks to Li-ion battery technology the power source has an operating time of typ. 9 hours, in continuous operation at maximum power typ. 4 hours. The mPack mini can be charged flexibly via the USB plug, e.g. on a laptop/desktop or via USB charger.



### NEW!

- ✓ 70 % smaller and lighter
- ✓ Operating time typ. 9 hours
- ✓ Charging time typ. 4 hours



Size comparison of HEINE mPack LL/mPack mini

### [ ORDERING INFORMATION ]

#### LoupeLight 2 Sets with HR Binocular Loupes

HR 2.5x optics with i-View on S-FRAME, LoupeLight 2, detachable yellow filter, mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, protective lenses, 2 sterilisable levers, retaining cord, cleaning fluid, cleaning cloth, carrying case

|                         |              |
|-------------------------|--------------|
| 340 mm working distance | C-008.32.450 |
| 420 mm working distance | C-008.32.451 |
| 520 mm working distance | C-008.32.452 |

#### MicroLight 2 Sets

MicroLight 2 on S-FRAME with mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, retaining cord and cleaning fluid in case

J-008.31.276

MicroLight 2 on Lightweight Headband with mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, in case

J-008.31.277

#### mPack mini

|   |              |
|---|--------------|
| mPack mini with Li-ion rechargeable battery and E4-USB plug-in transformer      | X-007.99.650 |
| mPack mini with Li-ion rechargeable battery, without E4-USB plug-in transformer | X-007.99.649 |
| Connecting cord SC 1 (1.6 m) for LoupeLight 2 or MicroLight 2 / mPack mini      | C-000.32.240 |

