






3M™ Littmann® Stethoscopes

Product Comparison Chart



	Electronic Stethoscope Model 3200	Electronic Stethoscope Model 3100	Master Cardiology™ Stethoscope	Cardiology S.T.C. (Soft Touch Chestpiece) Stethoscope	Cardiology III™ Stethoscope	Master Classic II™ Stethoscope	Classic II S.E. Stethoscope	Classic II Pediatric and Infant Stethoscope	Select Stethoscope	Lightweight II S.E. Stethoscope
										
Performance (1-10+ Scale)	10+	10+	10	9	9	8	7	7	6	6
Suggested Applications	Cardiology/ High Performance Bluetooth®	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Cardiology/ High Performance	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Physical Assessment, Diagnosis	Blood Pressure, Limited Physical Assessment	Blood Pressure, Limited Physical Assessment
Chestpiece	Durable Metal	Durable Metal	Polished Stainless Steel Single-sided	Matte Finished Stainless Steel, Single-sided with Soft Touch Insert	Machined Stainless Steel Two-sided Tunable	Plated and Polished Alloy Single-sided	Machined Stainless Steel Two-sided	Machined Stainless Steel Two-sided	Coated Alloy Single-sided	Metal/Resin Composite Two-sided
Headset Material	Aerospace Alloy	Aerospace Alloy	Aerospace Alloy	Aerospace Alloy	Aerospace Alloy	Aerospace Alloy	Aerospace Alloy	Aerospace Alloy	Aerospace Alloy	Aerospace Alloy
Length	27 in/69 cm	27 in/69 cm	22 in/56 cm 27 in/69 cm	27 in/69 cm	22 in/56 cm 27 in/69 cm	27 in/69 cm	28 in/71 cm	28 in/71 cm	28 in/71 cm	28 in/71 cm
Weight	185 grams	185 grams	185 grams	178 grams	175 grams	155 grams	125 grams	105 grams	115 grams	118 grams
Diaphragm Details, Diameter	51.0 mm	51.0 mm	Tunable 1.75 in/4.4 cm	Tunable, Teardrop 2.15 in x 1.9 in 5.4 cm x 4.8 cm	Two Tunable Diaphragms 1.7 in/4.3 cm	Tunable 1.75 in/4.4 cm	Tunable 1.75 in/4.4 cm	1.31 in/3.3 cm 0.75 in/1.9 cm	Tunable 1.75 in/4.4 cm	Tunable, Teardrop 2.15 in x 1.9 in 5.4 cm x 4.8 cm
Warranty	2 Years	2 Years	7 Years	6 Years	5 Years	3 Years	3 Years	3 Years	2 Years	2 Years