



Choose the stethoscope that helps you achieve your best.

3M™ Littmann® Stethoscopes offer superb acoustic performance, consistently high quality, and an outstanding value. Nearly all Littmann stethoscopes use tunable technology, a 3M invention that gives you the ability to hear different frequency sounds without turning over and repositioning the chestpiece.

Experience the power of tunable technology.



LIGHT CONTACT:
Hold the chestpiece with light pressure to hear low frequency sounds.



FIRM PRESSURE:
Press the chestpiece a little more firmly to hear higher frequency sounds.

Lightweight		Classic				Cardiology			Electronic				
6 Acoustic Performance Ideal for taking blood pressure readings and making limited patient assessments.		7-8 Acoustic Performance Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.				9-10 Acoustic Performance Designed to meet the needs of demanding clinicians. Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation.			10+ Acoustic Performance High-performance auscultation. Hear more of what you need to hear with amplification, ambient noise reduction, and optional connectivity.				
Lightweight II S.E. Stethoscope		Select Stethoscope		Classic II Pediatric and Infant Stethoscopes		Classic II S.E. Stethoscope	Classic III™ Stethoscope	Master Classic II™ Stethoscope	Cardiology III™ Stethoscope	Cardiology IV™ Stethoscope	Master Cardiology™ Stethoscope	Electronic Stethoscope Model 3100	Electronic Stethoscope Model 3200
Double sided chestpiece 118 grams Tunable diaphragm 2 year warranty Chestpiece: 2.12" / 5.4cm		Single sided chestpiece 120 grams Tunable diaphragm 2 year warranty Chestpiece: 1.75" / 4.4cm		Double sided chestpiece 105, 95 grams 3 year warranty Chestpiece: Pediatric 1.31" / 3.3cm, Infant 1.06" / 2.7cm, 0.75" / 1.9cm		Double sided chestpiece 135 grams Tunable diaphragm 3 year warranty Chestpiece: 1.75" / 4.4cm	Double sided chestpiece 150 grams Single-piece tunable diaphragm on each side of the chestpiece Convertible open bell 5 year warranty Chestpiece: Adult side 1.7" / 4.3cm, Pediatric side 1.3" / 3.3cm	Single sided chestpiece 160 grams Tunable diaphragm 3 year warranty Chestpiece: 1.75" / 4.4cm	Double sided chestpiece 175, 180 grams Available in 22" / 56cm and 27" / 69cm Two tunable diaphragms Dual-lumen tubing 5 year warranty Chestpiece: Adult side 1.7" / 4.3cm, Pediatric side 1.3" / 3.3cm	Double sided chestpiece 167, 177 grams Available in 22" / 56cm and 27" / 69cm Single-piece tunable diaphragm on each side of the chestpiece Dual-lumen tubing 7 year warranty Chestpiece: Adult side 1.7" / 4.3cm, Pediatric side 1.3" / 3.3cm	Single sided chestpiece 175, 185 grams Available in 22" / 56cm and 27" / 69cm Tunable diaphragm Dual-lumen tubing 7 year warranty Chestpiece: 2.0" / 5.1cm	185 grams 2 year warranty Chestpiece: 2" / 5.1cm	Bluetooth® wireless technology Recording and playback 185 grams 2 year warranty Chestpiece: 2" / 5.1cm

1-800-228-3957 | www.littmann.com

Bluetooth® is a registered trademark of Bluetooth SIG. 3M, Littmann, the L Littmann logo, L, Master Cardiology, Cardiology III, Cardiology IV, Master Classic II, Classic III, and the shape of the stethoscope chestpiece are trademarks of 3M. Used under license in Canada. Please recycle. Printed in U.S.A. © 3M 2016. All rights reserved. 70-2011-5462-5



