





Round Saw Blade TITANIUM

This newly developed saw blade is better than anything currently available. The main characteristic of this one of a kind product is the above average long life. It is nearly 400% in comparison with conventional saw blades. This is achieved not only due to high performance high speed steel (HSS) being used as the foundation material, but due to a special coating of Titanium Nitrate (TiN). This material combination makes it possible to cut very hard material, even metal. The TITANIUM Round Saw Blade fits all HEBU-Saws.




HB 8896-03
Round Saw
Blade TITANIUM
Ø 65 mm

HB 8895-03
Round Saw Blade
TITANIUM
Ø 50 mm



Non-stick Round Saw Blade

The technology and material of these saw blades correspond to the earlier HEBU blades proven thousands of times over. The difference, however, is that they feature a special non-stick coating, with the result that their anti-friction properties are improved at any speed ranges and the adherence of particles of dust is substantially reduced. These saw blades no longer stick and, as a result, make work easier and less tiring. They are especially suitable for all HEBU oscillating saws.



HB 8896-04
Non-stick
Round Saw Blade
Ø 65 mm

HB 8895-04
Non-stick
Round Saw Blade
Ø 50 mm



Round Saw Blade Ø 44 mm

For all HEBU Oscillating Saw Models

HB 8897-01 for plaster casts

HB 8897-02 for synthetic dressings, coated surface



Round Saw Blade Ø 50 mm

For all HEBU Oscillating Saw Models

HB 8895-01 for plaster casts

HB 8895-02 for synthetic dressings, coated surface



Round Saw Blade Ø 65 mm

For all HEBU Oscillating Saw Models

HB 8896-01 for plaster casts

HB 8896-02 for synthetic dressings, coated surface



Segment Saw Blade

Ø 65 mm

For all HEBU Oscillating Saw Models

HB 8892-01 for plaster casts

HB 8892-02 coated surface



Waisted Saw Blade

Ø 65 mm

For all HEBU Oscillating Saw Models

HB 8890-01 for plaster casts

HB 8890-02 coated surface



Deep Saw Blade

Ø 70 mm

For all HEBU Oscillating Saw Models

HB 8893-01 for plaster casts

HB 8893-02 coated surface

