



ASPEL

TREADMILL

B612

Speed	- speed range from 0 - 20km/h, settable by every 0,1km/h
Inclination range	- inclination range from 0 - 25%, settable by every 0,5%
Operating dimensions of band	- length - 1500mm - width - 500mm
ZERO START function	- function ZERO START available
Safety switch off button	- safety switch off buttons, stops treadmill. In standard are 2 pcs. on handrails and one external switch button for a medical personnel
In-built ECG module	- ECG signals: 12 standard leads in diagnostic mode - HR measuring range: 25 - 300bpm +/- 2% - HR resolution: 1bpm - monitoring range of ST section: +/- 3mm (10mm/mV) - ST wave change monitoring range: +/- 3mm, +/-0,2mm at 0,1mm resolution for 10mm/mV

ASPEL

TREADMILL

B612

In-built ECG module	<ul style="list-style-type: none">- ST wave change monitoring range: +/- 3mm, +/-0,2mm at 0,1mm resolution for 10mm/mV- sensitivity: 2,5/5/10/20mm/mV- record speed: 25/50/100mm/s Filters: <ul style="list-style-type: none">- network interferences: none, 50Hz, 60Hz- muscle interferences: none. 25Hz, 35Hz, 45Hz- isoline filter: none, 0,15Hz, 0,45Hz, 0,75Hz, 1,5Hz
Power supply	<ul style="list-style-type: none">- 230V, 50Hz
Max power consumption	<ul style="list-style-type: none">- 230V, 50Hz- below 10A
Max user weight	<ul style="list-style-type: none">- 180kg
Treadmill weight	<ul style="list-style-type: none">- 185kg
Dimensions	<ul style="list-style-type: none">- (LxWxH) 1999x747x1217mm
Safety norms	<ul style="list-style-type: none">- safety use: EN 60601-1, EN 957-1, EN 957-6, EN 60601-2-25- electromagnetic compatibility: EN 60601-1-2- medical device class: IIa (rule 10)- type of protection against electric shock (EN 60601-1): Class I- classification of the training equipment (EN 957-1): Class SA- application part (EN 60601-1): Type CF, resistant to defibrillation
Connection interface	<ul style="list-style-type: none">- RS-232 or USB - connection to PC which controls the work of a treadmill
Certificates	<ul style="list-style-type: none">- CE 0197- EN ISO 13485- MDD 93/42/EEC
Standard accessories	<ul style="list-style-type: none">- patient cable- disposable electrodes- abrasive paste