

NEW

Veinlite EMS PRO®

The world's bestselling vein access device just got better – now with integrated white exam light.



Veinlite®'s unique C-shaped design stretches the skin for easy needle insertion, prevents veins rolling and results in fewer IV access failures. Using Veinlite® saves time, avoids repeated needle sticks, decreases patient pain and anxiety and reduces wasted IV supplies.

- Clinically proven Veinlite® technology*
- NEW - Integrated exam light
- Robust design for use in the field
- Works on dark and light skin tones
- Uses any standard AA batteries
- Simple, single button operation

"They love it. They are constantly using it on all types of patients; pedis, geriatrics, and bariatrics. They also find that it helps locate veins on darker skinned patients. Thanks!"

Ruben Pena, Supply Chain Services Tech II, Emergency Department

SPECIFICATIONS

- 24 Bright Light Emitting Diodes (LEDs)
- Two LED colors optimized for vein imaging
- 12 Orange LEDs** for superficial veins
- 4 Red LEDs for deeper veins, darker skin
- 8 White LEDs for examining
- Light shield and pediatric adapter***
- Includes single use barrier covers
- Includes 2 AA batteries

TransLite, LLC Veinlite® Transilluminators
Designed and manufactured in the USA.



Patent Number US 7,874,698 B2 *Patent Number US 8,032,205 *Registered Trademark of TransLite LLC