

## SMART EMERGENCY BAGS

### POLYESTER



27153 27152 27150  
Shoulder straps included with medium and large bag. 4 pockets (5 for 27153) 1 or 2 on front, 2 on sides and 1 on back




27153

### POLYESTER WITH OUTER PVC COATING



27157 27158 27155  
Large pocket on top and shoulder straps included with medium bag



Safe metal buckles Large 10 mm zips Rubber bottom with feet Shoulder straps



Removable lining for easy cleaning

Large pocket on top

27157

Padded bottom and frontal pocket

**SMART BAGS: polyester or polyester PVC coated - empty**  
Professional emergency bags extremely capacious with internal adjustable compartments which permit personalized organization of the space. Provided with external pockets, internal divisions, and one (2 for code 27153) removable pouch with transparent window. All bags are made in strong, water resistant 600D polyester. Available also with outer PVC coating (27155-8). Provided with double yellow reflective strips. Internal lining is removable for easy cleaning. Waterproof rubber bottom to prevent water infiltration. May be carried on shoulder or hand.

Polyester	Polyester PVC coated	Model	Colour	Size	Internal modular spaces	External pockets
27150	27155	Small	red	45x28xh 28 cm	3	4
27151	-	Medium	red	55x35xh 32 cm	6	4
-	27157	Medium	red	55x35xh 38 cm	6	5
27152	-	Medium	blue	55x35xh 32 cm	6	4
-	27158	Medium	black	55x35xh 38 cm	6	5
27153	-	Large	red	65x35xh 35 cm	6	5

## EMERGENCY BAGS - POLYESTER OR POLYESTER PVC COATED

### POLYESTER

### POLYESTER WITH OUTER PVC COATING



27165 27167 27166



2 lateral pockets 22x18 cm and 11x18 cm

27167

Lateral pocket 32x18 cm

Front pocket 35x20 cm

- 27165 EMERGENCY BAG - polyester - red - 35x45xh 21 cm
  - 27166 EMERGENCY BAG - polyester PVC coated - red - 35x45xh 21 cm
  - 27167 EMERGENCY BAG - polyester PVC coated - blue - 35x45xh 21 cm
- Capacious bag, despite the compact size. Large central compartment, 1 front and 3 lateral pockets. May be carried on shoulder (braces) or hand (2 handles).

Made in strong, water resistant material (600D polyester). Available also with outer PVC coating (27166-7). Provided with yellow reflective strips and rubber bottom.