



The acoustic signal sounds when the defined platform capacity has been exceeded.



seca 862

Digital floor scale with cable remote control III

With its cable remote control (cable length: 2.3 m) for hand held, table or wall mounting, the battery-operated seca 862 is an extremely flexible, all-purpose scale. Energy-saving (up to 16,000 weighings per battery pack). In addition to its high capacity of 200 kg, the robust seca 862 has many other useful additional functions. For example, with the limit value function, it can be programmed to sound an alarm when a certain weight is exceeded. Among other things, this can be especially useful for pressure control in rehabilitation centres. The BMI function provides information on the nutritional condition of the patient.

- Capacity: 200 kg
- Graduation: 50 g < 150 kg > 100 g
- Functions: TARE, HOLD, BMI, limit value, automatic switch-off
- Optional: stand for cable remote display seca 472, RS232 interface seca 460, power adapter seca 447, switch-mode power adapter seca 400, carrying case seca 421

The raised rubber surface is robust and easy to clean.



seca 861

Digital floor scale with raised display III

With its large, easy-to-read, raised LCD display and the convenient tip-on and automatic switch-off, this scale is a particularly ideal solution for the medical sector. It is also very energy-saving (up to 16,000 weighings per battery pack) and robust: the casing can withstand heavy usage.

- Capacity: 200 kg
- Graduation: 50 g < 150 kg > 100 g
- Functions: tip-on automatic switch-on, automatic switch-off
- Optional: power adapter seca 447, switch-mode power adapter seca 400

With cushioned sides for absolutely safe transport.



Accessories

seca 421

Carrying case

With this carrying case of water-repellent material, the following seca scales can be safely transported: seca 862, seca 899 and seca 877.