

# BLOOD & INFUSION WARMER

## FT800

Effective, safe and easy to operate device for all medical departments:

- Blood transfusion / infusion during or after operation
- Routine blood transfusion and infusion in the clinical setting
- Enteral and parenteral nutrition
- Infusion to children and neonates
- Normal infusion in cold conditions

Standard infusion tube for medical fluids or blood

Microprocessor control

Actual heating temperature monitoring

One-touch-control

Holder for warming profile

Visual check of the infusion flow from the opening

Electric warming profile controlled by sensors

Cable to warming profile

Optimal heating up to the patient

# BLOOD & INFUSION WARMER

## FT800

### FEATURES

#### Heating up to the patient:

- IV tube is completely wrapped in, no heat loss

#### One-touch-control:

- Temperature automatically controlled
- Actual heating temperature in real time displayed

#### Safety features:

- Permanent running self-tests
- Double independent over-heating protections and automatic cut off
- Visual and acoustic alarm for high temperature / low temperature/sensor fault

#### Advanced dry heating technology:

- Quick warming up and effective heating

#### Open system:

- Accepts standard IV tube, no special disposables needed
- The most economical warming solution without extra consumable costs

#### Easy to operate:

- quick & simple installation
- 24 hours continuous operation

### TECHNICAL SPECIFICATIONS

Model:	FT800
Temperature setting:	39°C ex-works
Power supply:	a.c.100~240V/50~60Hz
Power consumption:	Max. 120VA
Type of protection against electric shock:	Class I
Degree of protection against electric shock:	BF Applied part; Defibrillation-protected
Degree of protection against ingress of liquids:	IPX2
Temperature accuracy:	±2.0°C
Overheat protection:	42°C/43°C
Low temperature alarm:	33°C
Warming up time:	From 20°C to 36°C approx.2 min.
Operating mode:	Continuous
Dimension (W*D*H):	85×65×175mm
Net Weight:	1.2kg
Warming profile:	length 1400mm; compatible to IV tube 3.5-5.0mm O.D. Contact the supplier for other options

\* Subject to technical modifications without notice



Infusion stand not included in delivery

**keewell**  
Medical Technology

Keewell Medical Technology Co.,Ltd.

Xiananyi Industrial Zone, Pingzhou, Nanhai, Foshan,  
Guangdong 528251, P.R.China

T.: +86 757 86260855 F.: +86 757 86260100  
marketing@keewell.cn www.keewell.cn

Distributor: