

Although rugged, the surface is easy to care for and can be cleaned with all commercially available disinfectants.




The measuring mat can be rolled up for space-saving storage.


## seca 210

Mobile measuring mat for babies and toddlers

Kind to the skin, washable and foldable: the light weight, space-saving solution for easy, precise measurement of the length of babies and toddlers while lying down. The fixed head piece and the sliding foot positioner make it simple to use.



- Measuring range: 10 – 99 cm
- Graduation: 5 mm
- With fitting for storing on wall



The scale's high-quality printing is especially tear-resistant. Measurements are precisely displayed and easy to read even after years of use.



For use with infants and toddlers, the measuring rod has head and foot pieces that slide.


## seca 207

Baby measuring rod with large callipers

For precise measuring when baby is lying down. The measuring rod can be firmly screwed to the changing table (wall mounting also possible) or used as a mobile rod. Large head and foot positioners make it easy to determine the length and the callipers fold for ease of storage. Made of aluminium for durability.

- Measuring range: 0 – 99 cm
- Graduation: 1 mm

A scale on each side permits measuring of head circumference and facial symmetry.



Precise and easy-to-read measurements made possible by clear red markings.

## seca 212

Measuring tape for head circumference of babies and toddlers

This measuring tape of non-stretch Teflon synthetic material measures both the head circumference (front side) and the symmetry of the face (reverse side). A practical dispenser (containing 15 tapes) ensures that they are stored safely.

- Measuring range: 5 – 59 cm
- Graduation: 1 mm