



Wall Transformer HEINE® EN50 UNPLUGGED

Intelligent charger for mPack UNPLUGGED



- ⚡ **Charger and wall-mount in one.** The wall transformer serves as an instrument holder and a charger at the same time.
- ⚡ **Intelligent charging system.** The instrument automatically switches off and goes into charging mode when placed in wall transformer. It is impossible to overcharge the mPack UNPLUGGED; the process ends automatically when charging is completed.
- ⚡ **Simultaneous charging.** The charging system charges the battery on the headband as well as an additional single spare battery.
- ⚡ **Compatible transformer (100 – 240V).** Multi-voltage power supply unit incl. 4 different mains adaptors.
- ⚡ **Registered Design.**

Wall transformer EN50 UNPLUGGED
(without mPack UNPLUGGED) 6V

X-095.17.310

⚡ Charger and wall-mount in one

⚡ Intelligent charging system

HEINE® mPack UNPLUGGED

Headband-mounted rechargeable battery

Rechargeable battery for OMEGA500 Indirect Binocular Ophthalmoscope UNPLUGGED and ML4 LED HeadLight UNPLUGGED (System voltage 6V).



- ⚡ **Lightweight with balanced weight distribution.** High level of comfort.
- ⚡ **Charge status indicator.** Display of remaining operating time.
- ⚡ Can be recharged by **wall- or plug-in transformer.**
- ⚡ **Direct operation possible while charging** through the use of the mobile plug-in transformer.
- ⚡ **Charging time:** typ. 2 hours.
- ⚡ **Operating time:** typ. 3.5 hours for ML4 LED HeadLight, typ. 2 hours for OMEGA500 UNPLUGGED, typ. 8 hours for OMEGA500 UNPLUGGED with LED.
- ⚡ **Intelligent charging system.** Automatic changeover from operating to charging mode (in the wall unit).
- ⚡ **Top safety standards.** Integrated charge and performance checks during operation. Lithium polymer technology. UL-compliant.
- ⚡ **Registered Design.**

⚡ Headband battery (rechargeable)

⚡ LED charge status indicator



[01]



[02]



mPack UNPLUGGED

mPack UNPLUGGED

X-007.99.665

Transformer UNPLUGGED for mPack UNPLUGGED [01]

X-095.16.316

Lithium-Polymer rechargeable battery for mPack UNPLUGGED

X-007.99.680

Extension Cord from plug-in transformer UNPLUGGED to mPack UNPLUGGED, 2 m [02]

X-000.99.668