

Ambu® ResCue Mask



Key Benefits

- The ResCue Mask is transparent to allow visual check of patient's condition
- Head strap for fixation and secure positioning. For use on both adults and children
- Reusable facemask that is easy to clean, which extends the product life
- Private labelling available
- Available in softcase

Ambu® ResCue Mask

The Ambu ResCue Mask is designed to provide protected rescue breathing for any patient (any age group). In CPR, mouth-to-mouth has been the traditional way of ventilation in an emergency situation. With the Ambu ResCue Mask first responders can provide immediate respiratory support without having any concerns about infectious diseases and the inherent unhygienic nature of the practice thus removing any natural reluctance to carry it out.

Features

The Ambu ResCue Mask helps protect rescuers in an emergency situation. It prevents direct contact with patient's mouth, nose and face and helps overcome hesitation to start resuscitation.

The mask is available with and without oxygen inlet. Thus an oxygen tube can be attached to the oxygen inlet for supplemental oxygen. The one-way filter of the mask can be removed in order to use the mask with a manual resuscitator.

The Ambu ResCue Mask comes in an original modern shell-shaped softcase.

Specifications

Product	
Size	
Dimensions	138 mm x 113 mm x 50 mm
Weight	Approximately 110g
Materials	
Mask	Polyvinylchloride
Softcase	Polyurethane, Ethylene Vinyl Acetate and Micro Fabric

Ordering specifications

Item no.	Ambu Rescue Mask
252 301 002	Ambu ResCue Mask – Adult mask with oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in black softcase
252 301 001	Ambu ResCue Mask – Adult mask with oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in white softcase
252 301 003	Ambu ResCue Mask – Adult mask with oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in blue softcase
252 301 006	Ambu ResCue Mask – Adult mask with oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in green softcase
252 301 004	Ambu ResCue Mask – Adult mask with oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase
252 301 102	Ambu ResCue Mask – Adult mask without oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in black softcase
252 301 101	Ambu ResCue Mask – Adult mask without oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in white softcase
252 301 103	Ambu ResCue Mask – Adult mask without oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in blue softcase
252 301 106	Ambu ResCue Maskv – Adult mask without oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in green softcase
252 301 104	Ambu ResCue Mask – Adult mask without oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase
Item no.	Ambu Rescue Mask Private Label
252 320 002	Ambu ResCue Mask – Adult mask with oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in black softcase
252 320 001	Ambu ResCue Mask – Adult mask with oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in white softcase
252 320 003	Ambu ResCue Mask – Adult mask with oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in blue softcase
252 320 006	Ambu ResCue Mask – Adult mask with oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in green softcase
252 320 004	Ambu ResCue Mask – Adult mask with oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase
252 320 102	Ambu ResCue Mask – Adult mask without oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in black softcase
252 320 101	Ambu ResCue Mask – Adult mask without oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in white softcase
252 320 103	Ambu ResCue Mask – Adult mask without oxygen inlet, elastic head strap and one-way valve with bacterial filter. Supplied in blue softcase
252 320 106	Ambu ResCue Mask – Adult mask without oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in green softcase
252 320 104	Ambu ResCue Mask – Adult mask without oxygen inlet, with elastic head strap and one-way valve with bacterial filter. Supplied in red softcase

*ResCue Mask with Oxygen Inlet is Rx Only in the USA

Logo
Organisation, company or sponsor logo. Maximum logo size: 70 x 40 mm.
Print
The softcase can be delivered in the most common Pantone colours.
The print colours can also be changed according to request.
Order
Minimum 1000 pieces.



The softcase for Ambu ResCue Mask comes in five standard colours.



The mask is packed inside the softcase with gloves and alcohol pad.

Ambu A/S
 Baltorpbakken 13
 DK - 2750 Ballerup
 Denmark
 Tel. +45 7225 2000
 Fax +45 7225 2053
 www.ambu.com