



Tongue-Blade Holder

For wooden tongue-blades



- **XHL Xenon Halogen technology.** Concentrated illumination of oral cavity and pharynx.
- **Light-guide for bulb.** Avoids dazzle.

Tongue-Blade Holder	2.5V XHL	3.5V XHL
for wooden tongue-blades, without handle	B-001.12.302	B-002.12.302
XHL Xenon Halogen spare bulb	X-001.88.037	X-002.88.049
Replacement light-guide for bulb [01]	B-000.12.308	

Tongue-Blade Holder

With blade ejector for use with HEINE disposable tongue-blades




- **XHL Xenon Halogen technology.** Concentrated illumination of oral cavity and pharynx.
- **Light-guide for bulb.** Avoids dazzle.
- **Blade ejector.** No need to touch used blades.
- **Only for HEINE disposable plastic tongue-blades.** Strong and hygienic.

Tongue-Blade Holder	2.5V XHL	3.5V XHL
with blade ejector, without handle	B-001.12.305	B-002.12.305
XHL Xenon Halogen spare bulb	X-001.88.037	X-002.88.049
Pack of 100 disposable plastic tongue-blades	B-000.12.304	

HEINE mini3000® Tongue-Blade Holder

With blade ejector for HEINE disposable tongue-blades



Can only be used with the mini3000 handle system.

- **Finger-tip pressure releases the used blade.** No need to touch the blade after use.
- **XHL Xenon Halogen technology.** Concentrated illumination of oral cavity and pharynx.
- **Only for HEINE disposable tongue-blades.** Strong and hygienic.

mini3000 Tongue-Blade Holder	2.5V XHL
Tongue-blade holder head with bulb, without handle [01]	D-001.74.100
complete with mini3000 battery handle* and 5 disposable blades, with batteries [02]	D-001.74.118
XHL Xenon Halogen spare bulb	X-001.88.037
Pack of 100 disposable plastic tongue-blades	B-000.12.304

* Available in black or blue. Please specify the colour you prefer when ordering. If not otherwise specified, we will supply black instruments.