

Ηλεκτρική Νοσοκομειακή κλίνη AC 505W

Μεταλλικός ανθεκτικός σκελετός με ξύλινες μετώπες, ξύλινα πλαινά και μεταλλικό σομιέ. Τροχήλατο με δυνατότητα ηλεκτρικής ανύψωσης της πλάτης και των ποδιών και καθ' ύψος ολόκληρης της κλίνης. Επίσης διαθέτει επαναφορτιζόμενη μπαταρία για χρήση του μηχανισμού σε περίπτωση διακοπής ρεύματος.



Fig. 5: Angles of the bed base

Bed base:	94 x 200 cm
Lowest position of bed base	40 cm
Highest position of bed base	80 cm
Overall dimensions	105 x 214 cm
Idle weight	Bed: 84 kg Lifting pole: 6,4 kg Side rails: 9 kg Accessories: 0,25 kg
Maximum permitted load:	<u>Bed</u> Max. weight of the user: 135 kg mattress 20 kg accessories 15 kg Total 170 kg Lifting pole and triangle grip: 80 kg
Maximal load of the lifting pole:	80 kg
Maximum load total	170 kg
Heaviest sections of the bed: Head/foot of the bed	Weight : 20 kg
Back section :	70° (see figure 5)
Foot section :	0-30° thigh section (see figure 5) 0-20° leg section : 5 notched stages (see figure 5)
Trendelenburg / Reverse Trendelenburg	± 11°
Electrical category	Class II
Noise level	Less than 65 dB
Maximum permitted mattress height:	12 to 18 cm (Wooden side rails) 12 to 19 cm (Metal side rails)
Mattress dimension	90 x 200 cm