

iS3 Video Laryngoscope


Economical Tool,
Easier Intubation Experience



Innovation for Simplifying Intubation

At Insighters, innovation for clinical demands is a core driving force in our R&D process.

The iS3 Video Laryngoscope, as a new member of Insighters i-Scopes product family, its inherited innovative technologies and structure designs that have been widely accepted by world-wide users.



"Insight Slider" design

Support physicians to perform intubation for all sizes of patients by just one handle, Easy&Convenient



Real-time Anti-fog

Instant and real time anti-fog function, no need coating or waiting for preheating



Multi-media & Inside Menu

Integrated photography and video recording, and built-in guide video, in-door/out-door screen mode, easy to use and learn.

Disposable & Reusable Blade

Blade L
(Large Adult)



Blade M
(Adult)



Blade S
(Child)



Blade SS
(Neonate)



PRODUCT SPECIFICATIONS

iS3 Host

Display	3.5" LCD display, touch screen
Battery	Rechargeable lithium battery
Working Time	240 minutes

iS3 Handle


Size	120 mm x 150 mm x 36 mm
Weight	80g
Camera	CMOS
Light source	White LED

iS3 Blades

Material	Polycarbonate
ET Size	≥ 2.5 mm

Portable, Powerful, Economical

- Light and easy to hold with perfect center of gravity when do intubation.
- Easy and convenient to carry where needed.
- Screen brightness adjustable for different location.
- Host and handle detachable, easy cleaning and disinfection.
- Rechargeable lithium battery, economical and durable.
- Support Macintosh 3 & 4 blades and insighters blades, for routine and difficult airway.



Shenzhen Insighters Medical Technology Co., Ltd
Floor 13, HTM Building, NO.3 Ganli Road, Longgang
District, Shenzhen, Guangdong, China.
Tel: +86-755-28218609 ext.861
Email: intl-market@insighters.cn
www.insighters.cn

