



seca 634
EMR ready scale for heavy persons with large step on surface, cable-link remote control and high load-bearing capacity.

Page 54



seca 869
Robust flat scale with cable remote display, BMI, HOLD and TARE for mobile use.

Page 55



seca 874
Electronic flat scale for mobile use with push buttons and double display.

Page 56



seca 876
Robust, very light flat scale with integrated display, tip-on automatic switch-on and mother/child function for mobile use.

Page 58



seca 874 dr
Doctor's scale with customizable label.

Page 58



seca 813
Electronic flat scale with very high bearing capacity of up to 200 kilograms.

Page 60



seca 803
Electronic flat scale with very high bearing capacity of up to 200 kilograms.

Page 60



seca 762
Mechanical flat scale with large dial and fine 500 g graduation.

Page 61



seca 750
Mechanical flat scale with robust, powder-coated steel casing. Kg/lbs and lbs version available.

Page 61

TECHNICAL DATA

| | | | | | | | | | |
|--------------------------------|---|---|--|---|--|---|---|---|--|
| Capacity | 360 kg / 800 lbs | 250 kg / 550 lbs | 200 kg / 440 lbs | 250 kg / 550 lbs | 200 kg / 440 lbs | 200 kg / 440 lbs | 150 kg / 330 lbs | 150 kg, 320 lbs, 150 kg / 320 lbs | 150 kg, 320 lbs, 150 kg / 320 lbs |
| Graduation | 50 g / 0.1 lbs | 100 g < 150 kg > 200 g 0.2 lbs < 330 lbs > 0.5 lbs | 50 g < 150 kg > 100 g 0.1 lbs < 330 lbs > 0.2 lbs | 100 g < 150 kg > 200 g 0.2 lbs < 330 lbs > 0.5 lbs | 50 g < 150 kg > 100 g 0.1 lbs < 330 lbs > 0.2 lbs | 100 g / 0.2 lbs | 100 g / 0.2 lbs | 500 g, 1 lbs, 500 g / 1 lbs | 1 kg, 1 lbs, 1 kg / 1 lbs |
| Dimensions (WxHxD) | 640 x 53 x 600 mm 25.2 x 2.1 x 23.6" | 321 x 60 x 356 mm 12.6 x 2.4 x 14" | 321 x 60 x 361 mm 12.6 x 2.4 x 14.2" | 321 x 60 x 356 mm 12.6 x 2.4 x 14" | 321 x 60 x 361 mm 12.6 x 2.4 x 14.2" | 433 x 47 x 373 mm 17 x 1.9 x 14.7" | 316 x 37 x 326 mm 12.4 x 1.5 x 12.8" | 305 x 120 x 470 mm 12 x 4.7 x 18.5" | 317 x 117 x 477 mm 12.5 x 4.6 x 18.8" |
| Dimensions/platform (WxHxD) | 560 x 55 x 560 mm 22 x 2.2 x 22" | 280 x 60 x 278 mm 11 x 2.4 x 10.9" | 288 x 60 x 280 mm 11.3 x 2.4 x 11" | 288 x 60 x 278 mm 11.3 x 2.4 x 10.9" | 288 x 60 x 280 mm 11.3 x 2.4 x 11" | 400 x 45 x 300 mm 15.7 x 1.8 x 11.8" | 295 x 35 x 305 mm 11.6 x 1.4 x 12" | 260 x 70 x 265 mm 10.2 x 2.8 x 10.4" | 290 x 70 x 290 mm 11.4 x 2.8 x 11.4" |
| Cable length | 2.5 m / 98.4" | 2 m / 78.7" | | | | | | | |
| Weight | 14.4 kg / 31.7 lbs | 4.5 kg / 9.9 lbs | 4.1 kg / 9 lbs | 4.2 kg / 9.3 lbs | 4.1 kg / 9 lbs | 2.9 kg / 6.3 lbs | 1.5 kg / 3.3 lbs | 3.5 kg / 7.7 lbs | 3.5 kg / 7.7 lbs |
| Power supply | Batteries, power adapter optional | Batteries, power adapter optional | Batteries | Batteries | Batteries | Batteries | Batteries | | |
| FUNCTIONS | | | | | | | | | |
| TARE | • | • | | | | | | | |
| Pre-TARE | • | | | | | | | | |
| Mother/child function | | | • | • | • | | | | |
| HOLD | • | • | | | | | | | |
| Auto-HOLD | • | | • | | • | • | • | | |
| Step-off | | | | | | | | | |
| BMI | • | • | | | | | | | |
| kg/lbs switch-over | • | • | • | • | • | • | • | | |
| Auto-CLEAR | • | | | | | | | | |
| Damping | • | | | | | | | | |
| Reset-to-zero function | | | | | | | | • | • |
| Tip-on automatic switch-on | | | | • | | • | • | | |
| RESET | • | | | | | | | | |
| Automatic switch-off | • | • | • | • | • | • | • | | |
| SEND/Auto-SEND | • | | | | | | | | |
| PRINT/Auto-PRINT | • | | | | | | | | |
| OPTIONAL | | | | | | | | | |
| Interface module | seca 452, seca 456 | | | | | | | | |
| Measuring rod | | seca 217 | seca 217 | seca 217 | | | | | |
| Stand for cable remote display | seca 472 | | | | | | | | |
| Adapter element | | seca 437 | seca 437 | seca 437 | | | | | |
| Power adapter | seca 447, seca 400 | seca 447, seca 400 | | | | seca 447, seca 400 | | | |
| PC interface cable | | | | | | | | | |
| Carrying case | | seca 409, seca 414, seca 421 | seca 409, seca 414, seca 415 | seca 409, seca 414, seca 415 | | | seca 415 | | |
| Individualizable label | | | | | seca 478 | | | | |
| AREAS OF APPLICATION | | | | | | | | | |
| Mobile weighing | | • | • | • | | | | | |
| EMR READY | • | | | | | | | | |

