



HEINE rechargeable batteries for BETA® Handles



Rechargeable batteries	3.5V NiMH	3.5V LI ION	3.5V LI ION L
for BETA rechargeable handles	X-002.99.382	X-007.99.381	X-007.99.383
	[01]	[02]	[03]

Please note: The Li-ion L battery is not compatible with BETA R, TR and NT 3.5V handles and bottom inserts. For information on HEINE Laryngoscope Handles, see pages 085 – 089.

- All full-size BETA handles can be updated to BETA L 3.5V (except BETA SLIM) handles by substituting the Li-ion L 3.5V battery and the BETA L bottom insert.
- Max. 4 hours recharge time in NT300 table charger with NiMH. 2 hours with Li-ion.
- Retro-compatible. All rechargeable batteries can also be charged in existing NT200 and Nicatron chargers.



BETA Belt Clip	X-000.99.008
-----------------------	---------------------

HEINE Large Battery Handle



Dimensions: 165 mm x 37 mm dia. Weight: 212g without batteries, 505g with batteries.
A large handle for 2 alkali-manganese batteries IEC LR20/D cells.
Designed for countries where larger batteries are still preferred.
(Only suitable for certain instrument cases.) Batteries not included.

Large battery handle	2.5V
for 2 alkali-manganese batteries, size D	X-001.99.120

HEINE BETA® SLIM Battery/Rechargeable Handle



Dimensions: 135 mm x 19mm dia. Weight: 180g.
A compact handle with automatic connector for those who prefer a slim handle.
Rechargeable in NT300 charger fitted with the small adaptor ring.
Batteries not included.

BETA SLIM battery / rechargeable handle	2.5V	3.5V NiMH
BETA SLIM Battery Handle	X-001.99.105	
BETA SLIM Rechargeable Handle with 3.5V M3Z NiMH battery		X-002.99.471
Rechargeable battery 3.5V M3Z NiMH [01]		X-002.99.106

[01]