



**Laerdal**<sup>®</sup>

helping save lives



## Fundamental Pediatric Training

### Nursing Kid

Nursing Kid realistically represents a six year old and is ideal for clinical training of the care and management of pediatric patients. Nursing Kid paired with the SimPad<sup>®</sup> System helps facilitate delivery of consistent simulations into nursing curricula.



[laerdal.com](http://laerdal.com)

Visit us to learn more about  
Nursing Kid and SimPad<sup>®</sup>



# Nursing Kid

## PATIENT CARE

- Irrigation of the eye and ear (simulated)
- Nasal packing
- Ear, eye, nose and mouth care
- Gastric lavage and gavage
- Full range of motion for realistic patient handling
- Dressing and bandaging techniques possible
- Pediatric transfer techniques possible

## AIRWAY

- Insertion, securing and care of endotracheal tubes
- Insertion and suctioning of oropharyngeal and nasopharyngeal airways
- Tracheostomy care and tracheal suctioning
- NG tube insertion, care and removal
- Various oxygen delivery procedures

## PULSE

- Manually generated carotid pulse

## INTERCHANGEABLE GENITALIA

- Interchangeable male and female genitalia with connector valves and urinary and colon reservoirs

- Complete urinary catheterization
- Will retain indwelling or straight catheter
- Enema procedures may be performed

## IV/IM TRAINING

- Articulating IV arm with replaceable skin and infusible vein system allows peripheral intravenous therapy and site care
- Venipuncture possible in the antecubital fossa and dorsum of the hand
- Accessible veins include median, basilic and cephalic
- Thigh (bilateral), deltoid, and gluteal intramuscular injections (unilateral)

## SOUNDS - SIMPAD® SYSTEM

- Heart sounds synchronized with ECG
- Auscultated lung sounds synchronized with breathing, 0 - 60 BPM
- Individual lung sound selection
- Normal or abnormal bowel sounds
- Vocal sounds: Computer-generated sounds, recorded vocal sounds and real-time voice input via headset
- User generated vocal sounds

## SIMPAD® SYSTEM FUNCTION

- Handheld, intuitive touchscreen remote for easy 'pick up and play' experience
- Mobile - teach anywhere
- Operate on-the-fly or utilize scenarios and Themes for consistent simulation training
- Time stamped activities, vital signs, and instructor comments are captured in the event log
- The SimPad log files can be viewed on the device or on a PC using Session Viewer for post-simulation reflection and debriefing
- Upload self-authored scenarios and Themes, or download pre-programmed scenarios directly from SimStore

## PATIENT MONITOR (OPTIONAL)

- Touchscreen simulated patient monitor provides concise clinical feedback for physiological parameters.
- The monitor's color screen is configurable and provides multiple simulated parameters, each presenting multi-level alarms.
- Simulated parameters include HR, ECG, SpO2, BP, RR, Temperature, and etCO2

## ORDERING INFORMATION

### 350-05050 Nursing Kid Light Skin

Includes: Full Body (six year old) SimPad-Capable Manikin, Airway Lubricant, Simulated Blood, Hospital Gown, and Carry Case.

### 277-00001 Pediatric Wound Modules



### 200-30001 SimPad® System

Includes: SimPad, Link Box, 2 AC Power Supplies, Rechargeable Lithium-Ion Battery, Wrist Strap for SimPad, Manikin Strap for Link Box, Manikin Adapter Cable, SimPad Sleeve, USB Cable, Ethernet Cable and Directions for Use.

### 400-09201 Tablet Simulated Patient Monitor (11.5")

### 400-29301 All in One Simulated Patient Monitor (23" panel)

## Operation

Nursing Kid is controlled by the SimPad® System and can be utilized with an optional simulated patient monitor.



\*\* SimPad Unit required for operation. Sold separately \*\*

## SimCenter

Visit [www.laerdal.com/simcenter](http://www.laerdal.com/simcenter) to learn how to optimize your simulation training using SimCenter™ and Nursing Kid



SimStore™

Access to validated, standardized content from industry experts

SimManager™

A suite of online tools that allows you to fully optimize your entire simulation center operation

SimView™

A debriefing solution giving you the ability to study simulation events more closely

## Services

Laerdal offers comprehensive Educational Services to support your curriculum, and Technical Services to meet your product servicing needs. Whether you are setting up a one-room simulation site or a multi-bed simulation center, Laerdal can help.

Visit [Laerdal.com/us/NursingKid](http://Laerdal.com/us/NursingKid) for ordering information and a complete list of product accessories.

