

NObreath® comparison chart

		
Adult breath test mode	✓	✓
Child breath test mode	✓	✓
Ambient breath test mode	✓	✓
Demo breath test mode	✗	✓
Verification of a successful test	✗	✓
Verification of an unsuccessful test	✗	✓
Ability to store patient profiles	✗	✓
Ability to show results on a graph	✗	✓
Bluetooth connectivity	✗	✓
PC connectivity	✗	✓
Internal mass flow sensor	✗	✓
Breath test animations for motivation	✗	✓
0-500ppb measurement range	✗	✓
3l/m breath test flow rate	✓	✓
5 years sensor operating life	✓	✓
SteriTouch® technology	✗	✓
Rechargeable battery	✗	✓
Charging/syncing dock	✗	✓
Plug & play serviceable parts	✗	✓
Pre-calibrated/changeable sensor	✗	✓
NObreath® forum	✓	✓
Free FeNOchart™ patient management software	✗	✓

Bedfont® Scientific Ltd.

Station Road, Harrietsham, Maidstone, Kent, ME17 1JA, England

Tel: +44 (0)1622 851122, Fax: +44 (0)1622 854860

Email: ask@bedfont.com www.bedfont.com

© Bedfont® Scientific Limited 2017

Issue 1 – September 2018, Part No: MKT507

Bedfont® Scientific Ltd. reserves the right to change or update this literature without prior notice. Registered in England and Wales. Registered No: 1289798