



Throat Scope®



Manufactured by: ThroatScope Pty Ltd. 1/48 Kalaroo Road, Redhead NSW 2290 Australia

www.throatscope.com

Throat Scope

Illuminated tongue depressor single use blades.

Features

- Light and tongue depressor in one which can simultaneously depress the tongue and illuminate the mouth
- Slim, lightweight unit
- Light automatically activated
- One handed operation
- Reusable light
- Single use disposable blade
- Powerful long life 3 x LED lights

Operational Instructions

Thank you for purchasing a 50 or 100 pack Illuminated Tongue Depressor Single Use Blades consumable pack. Please read the following operational instructions carefully before use.

Note: Illuminated tongue depressor light sold separately.

Precautions and Warnings

- Disposable blade is single use. **DO NOT REUSE.** Do not attempt to clean the blade or place in the dishwasher. Please dispose of blade after use.
- Between use, clean the light by wiping with 70% isopropyl alcohol to reduce the risk of cross-contamination.
- For oral use only. The device is intended to depress the tongue and illuminate the mouth.
- The device is intended to be used by healthcare professionals. In the event that you are uncertain or uncomfortable about using the device, see a healthcare professional.
- Do not bend or bite the light.
- Do not attempt to disassemble the light.
- Dispose of Light once battery is flat.
- Keep out of the reach of unsupervised children.

Directions

1. Insert the single use blade into the light. The light will automatically activate when inserted.
2. Insert the device into the oral cavity to use.
3. When finished, remove the blade by pushing on the blade grips and dispose of the blade. The light will automatically turn off when removed.
4. After use, clean the light by wiping the entire light surface thoroughly with a 70% isopropyl alcohol wipe.

