



Armchairs • Accessoires



TECNOLOGIE DI PESATURA



# MODENA



## Movements:

- Backrest lifting/lowering
- Leg part lifting/lowering
- Simultaneous backrest/leg part lifting/lowering
- Variable height lifting /lowering allows the patient to get out the armchair



## Features:

**BASE:** new ergonomic design Ergonomic design of new conception that ensures greater stability to the armchair with less obstacles; the particular design, fluid and dynamic, makes it quicker and easier cleaning and sanitizing

**ARMREST:** Padded, shaped, manually adjustable, opening facing the outside to allow the patient to sit easily

**PADDING:** Equipped with two padded headrest, adjustable and removable. The interchangeable mattress is made with differentiated density padding.

**FOOTREST:** With manual adjustment thought to prevent cramps (available electric version).

**WHEELS:** with breaking system

Three section seat, each adjusted independently, operated by push-button panel. Layered padding to give the patient greater comfort; washable, fireproof, atoxic, removable, can be disinfected, eco-leather covering; with interchangeable mattress and two removable headrests. Arms with integrated push-button panel holder and synchronized, ergonomic movement keeping the elbow in the correct position based on the different backrest inclinations.

## DATI TECNICI

Size support plane	210 x 56 cm
Size with armrests	210 x 80 cm
Minimum Height	55 cm
Maximum Height	90 cm
Tilt Seat	0-13°
Backrest Angle (with respect to the seat tilt)	0-80°
Leg-rest Angle (with respect to the Sitting)	0-20°
Trendelenburg	13°
Footrest excursion	23 cm
Wheels diameter	12,5 cm
Safe Working Load	200 Kg
CND Code	Z129008
Progressive Code	1210547



**Pushbutton:** the armchair is equipped with a pushbutton panel both for the patient and the operators in order to reach each movement of the various sections; can be located through a hook when not use. It is protected against any liquid infiltrations by a membrane.

## OPTIONAL ACCESSOIRES

ACCESSOIRES		MODENA	IDEA
302570	PUSHBUTTON TRAY	✓	✓
327000	SERVING TROLLEY IN LAMINATED HEIGHT ADJUSTABLE	✓	✓
327600	SERVING TROLLEY IN STAIN-PROOF ABS HEIGHT ADJUSTABLE	✓	✓
310580/59	OBJECT-HOLDER TABLE	✓	
310580	OBJECT-HOLDER TABLE WITH IV POLE	✓	
346007	IV POLE	✓	
ATT/FL	IV POLE ATTACHMENT	✓	✓
PRT	ROLL STAND	✓	•
BR/AV	HEIGHT ADJUSTABLE ARMREST	✓	
C/LOM	LUMBAR PILLOW	✓	
PED/EL	ELECTRIC FOOTREST	✓	•
LE/SB	LEVER FOR FAST UNLOCKING	✓	
AV/P	VARIABLE HEIGHT FOR ARMCHAIR	✓	•
311100	COURTESY LAMP		✓
311110	FLEXIBLE PUSHBUTTON TRAY		✓
311120	ADDITIONAL PADDING SUPER COMFORTABLE		✓
311140	ADJUSTABLE IV POLE WITH TRAY		✓
311142	IV POLE WITHOUT TRAY		✓

• STANDARD



### COLOURS AVAILABLE



## RELAX



FLY



DAY



EURO



LUCCA



MADRID



NAPOLI



ROSE

5.9 s.r.l. was established in 2014 from the WAGER, CHALLENGE and COMMITMENT of professionals who wanted to make the knowledge gained during a work experience at Tassinari Balance (worldwide known brand in the weighing sector) available to patients and users, through its production of bed scales and armchair scales for medical application.



5.9 s.r.l. uses qualified professionals to manage their Customer / After Sales Service with excellent skills. Expertise and know-how are the main features of the service, as they are an integral part of the company's commitment that dedicates the same attention to maintenance and repairs as it does to production.

The experienced technical staff is able to maintain the efficiency of each instrument offering a highly professional service, using quality spare parts. This service is available widely throughout the national territory (including islands), also through authorised Technical Assistance Centres (CAT) and the HELP-DESK service is always in operation for any requirements

Cod. CAT/POLTR Rev.01 26.10.16  
"L'azienda si riserva di apportare ai prodotti, senza alcun preavviso, tutte le modifiche tecniche e/o estetiche ritenute opportune".



## 5.9 s.r.l. CARE WEIGHTING SYSTEM

Via Ischia, 2 - 40017  
San Matteo della Decima (Bo) - Italy

T. +39 051 0546950  
F. +39 051 19901816  
info@5dot9.it