



PARI COMPACT

Quality at an affordable price

- Nebuliser system for diseases of the lower airways, such as COPD, Asthma or Bronchitis
- Proven efficacy: includes **PARI LC PLUS**[®] nebuliser with mouthpiece and mask

Technical data:

Operating pressure:	1.0 bar
Weight:	1.2 kg
Dimensions (w x h x d)	10 x 14 x 16 [cm]
Warranty:	2 Years

Aerosol Performance:

TOR [mg/min]	418
MMD [μ m]	3.9
Mass percentage < 5 μ m	64%

Measured with the Malvern Mastersizer X (according to Fraunhofer-Model) at 21 °C and 50% relative humidity. Nebulized medium: 0.9% NaCl (5 ml). Inspiratory flow 20 l/min

PARI COMPACT and its accessories



PARI COMPACT

Nebuliser system 230-240 V~/50Hz with PARI LC PLUS® Nebuliser, Child Mask, Mouthpiece and Connection Tubing

Item No. 052G1004



PARI LC PLUS® Nebuliser

Annual Replacement Pack. Consists of LC PLUS® Nebuliser, Connection Tubing (2.1 m) and Filter Holder with Filter.

Item No. 022G8107



Filter Holder with Filter

Item No. 041B6500



PARI Child Mask (Soft)

with Elastic Band

Item No. 041G0741



PARI Adult Mask (Soft)

with Elastic Band

Item No. 041G0740



Connection Tubing (f/f)

Item No. 041B4592



PARI Medical Ltd
The Old Sorting Office
Rosemount Avenue
West Byfleet, Surrey KT14 6LB
Tel: 01932 341122
Fax: 01932 341134
email: infouk@pari.eu

www.parimedical.co.uk

Specialists in effective inhalation

