



BLOOD PRESSURE MANAGEMENT

Full portfolio of blood pressure devices and cuffs that allow true standardisation and help achieve excellence in BP management.



Sphygmomanometers

BLOOD PRESSURE MANAGEMENT

ProBP® 2000 Digital Blood Pressure Device

A fast, accurate and cost-effective automated blood pressure device for healthcare professionals looking for both clinical and economical efficiency.

- > Combines two of our innovative technologies—Welch Allyn SureBP technology and Welch Allyn FlexiPort® Blood Pressure Cuff technology—to enhance blood pressure tests for caregivers and patients
- > Blood pressure readings are captured more comfortably on inflation and take about 20 seconds
- > Compatible with Reusable FlexiPort cuff sizes, ranging from child up to large adult and thigh (15-55cm range)
- > One-button touch triggers the device to start capturing on cuff inflation
- > Digital LCD screen displays BP readings and pulse rate
- > Works with AA batteries or with optional AC adaptor

2000 ProBP 2000 with adult Reusable FlexiPort BP cuff (size 11) and 4 AA batteries



7670-01

767 Wall and Mobile Aneroids

The versatile 767 Wall and Mobile Aneroid, with improved ergonomics and a new modern design, is the smart choice for accurate and reliable readings and superior lifetime value.

- > Lifetime calibration warranty
- > Laser-engraved dial face provides unmatched accuracy
- > Recessed dial face and high contrast red pointer increase visibility and reduce parallax errors
- > New 5-leg stand provides improved stability and manoeuvrability
- > Latex and mercury-free

- 7670-01** 767 Wall Aneroid with Reusable Adult FlexiPort Cuff
- 7670-10** 767 Mobile Aneroid with Reusable Adult FlexiPort Cuff
- 7670-16** 767 Desktop Aneroid with Reusable Adult FlexiPort Cuff



7670-16

7670-10