

## RS01 Sleep Apnea Screen Meter

The RS01 device is suitable for the one who suffers with such diseases as SAHS (sleep apnea-hypopnea syndrome), Its functions including patient's sleep apnea monitoring and SpO<sub>2</sub> monitoring, which can be applied in hospital or at home.

### Features

- ◆ SpO<sub>2</sub> and PR value display, pulse sound and bar graph indication.
- ◆ Pulse waveform and nose flow waveform display.
- ◆ Alarms for low power, finger out and exceeding limits.
- ◆ Battery status indication.
- ◆ Real-time clock.
- ◆ Timing power on/off.
- ◆ Storage of multiple cases.
- ◆ Continuous 12 hours data record.
- ◆ Uploading data via USB port.
- ◆ PC analysis software.







RS01

## RESTING ECG

## PM10 Portable ECG Monitor

### Features

- ◆ 1.77" color TFT-LCD.
- ◆ Quickly to check ECG once finger touching, convenient to operate.
- ◆ Accurate conclusion can be obtained immediately after measuring.
- ◆ The built-in rechargeable lithium battery can continuously record for up to 500 cases in full charging.
- ◆ Bluetooth transmission.
- ◆ ECG data can be saved to CLOUD Platform permanently, convenient to check and analyze.
- ◆ Historical data can be printed at any time, convenient to diagnose.



PM10