

HOLTER



TLC5000/TLC9803/TLC6000 Dynamic ECG Systems

Recorder features

- ◆ Strong in anti-interference and anti-vibration capability.
- ◆ Real-time clock(only for TLC5000, TLC9803 and TLC6000).
- ◆ With the functions of waveform preview, record review and event marking(only for TLC5000 and TLC9803).
- ◆ Record ECG waveform in detail by high-accuracy and high-sampling frequency recorder.

Recorder parameters

	TLC5000	TLC9803	TLC6000
Lead	12-lead	3-lead	12-lead
Memory	MicroSD(1GB)	MicroSD(1GB)	MicroSD(1GB)
Record time	24h	24h	24-48h
Power supply	Two "AA" batteries	Two "AA" batteries	One "AAA" battery
Screen	1.77" TFT	1.77" TFT	0.96" OLED
Interface	USB 2.0	USB 2.0	USB 2.0
Dimension	111mm(L)×60mm(W) ×25mm(H)	111mm(L)× 60mm(W) ×25mm(H)	80.5mm(L)×59.5mm(W) ×22mm(H)



TLC9803/TLC5000



TLC6000