

<b>78812 (EU) / 78814 (UK)</b>	KleenSpec Corded Illumination System—Light Handle with Integrated On/Off Switch and Transformer
<b>78800</b>	Light Handle with On/Off Switch Only

## KleenSpec Corded Illumination System: Accessories

<b>73412 (EU) / 73414 (UK)</b>	Transformer only for 580 and 586 Series Illumination System, 1.5 m cord
<b>08800-U</b>	Replacement 6 V Halogen Lamp for No. 78812 / 78814

## KleenSpec® Vaginal Specula

Welch Allyn's single-use specula eliminate the time and expense of reprocessing and reduce the chance of cross-contamination. Each speculum is designed and manufactured to deliver durability and reliability.

- Newly designed light pipe provides a homogeneous spot and eliminates glare-back
- New wider handle provides improved ergonomics, better balance, and is much easier to manipulate during exams
- Smooth, molded acrylic will not pinch or bind tissue for greater patient comfort during the exam
- The new 590 Series KleenSpec Vaginal Speculum is compatible with the new cordless illumination system (79912/79914) as well as the illumination system (78812/78814)
- Pederson Style Corded

<b>59000</b>	KleenSpec 590 Series Disposable Vaginal Specula, Small (24/box; 4 boxes/case)
<b>59001</b>	KleenSpec 590 Series Disposable Vaginal Specula, Medium (24/box; 4 boxes/case)
<b>59004</b>	KleenSpec 590 Series Disposable Vaginal Specula, Large (18/box; 4 boxes/case)



**59000**  
590 Series Disposable Vaginal Specula, Small



**59001**  
590 Series Disposable Vaginal Specula, Medium



**59004**  
590 Series Disposable Vaginal Specula, Large