



WOMEN'S HEALTH

Single-use vaginal specula and cordless illumination system designed to help women's health providers improve patient comfort and safety.



Vaginal Speculas and Illumination System


WOMEN'S HEALTH

KleenSpec® 800 Cordless Illumination System

Designed with enhanced unique new features, including bright, white light with our new LED light sources.

- > Advanced LED light source provides improved visualisation of the examination area by supplying brighter, whiter light
- > No cord to get in the way during procedures
- > Improved durability to better withstand cleaning & disinfecting

- 80015** Complete KleenSpec 800 Disposable Vaginal Specula Cordless Illuminator with Charging Station
- 80000** KleenSpec Disposable Vaginal Specula Cordless Illuminator
- 59010** Disposable Sheath for the KleenSpec Cordless Illumination System



KleenSpec Vaginal Specula—Premium 590 Series

Welch Allyn's single-use specula are designed to eliminate the time and expense of reprocessing and can help reduce the chance of cross-contamination.

- > Now in every size needed for all types of patients—safe examination with virginal/paediatric patients as well as post-menopausal and post-hysterectomy women
- > Look for the distinctive blue, green, white and yellow coloured handles to easily identify sizes
- > Newly designed light pipe provides a homogeneous spot and eliminates glare
- > New wider handle provides improved ergonomics, better balance, and easy manipulation during examinations
- > Smooth, moulded acrylic will not pinch or bind tissue for greater patient comfort
- > Compatible with both KleenSpec illumination systems (cordless and corded)
- > Industry's first extra small size disposable vaginal specula

- 590XS** KleenSpec 590 Disposable Vaginal Specula, Extra Small (case of 96)
- 59000** KleenSpec 590 Disposable Vaginal Specula, Small (case of 96)
- 59001** KleenSpec 590 Disposable Vaginal Specula, Medium (case of 96)
- 59004** KleenSpec 590 Disposable Vaginal Specula, Large (case of 72)

KleenSpec Vaginal Specula—Premium 590 Series with Smoke Tube

For smoke evacuation during electrosurgical procedures.

- > Built-in smoke evacuation channel to assure a clear view
- > Advanced LED light source eliminates positioning of an external light source
- > Convenient single-use disposable for patient safety
- > Excellent visibility—transparent blades allow full visualisation of the examination area
- > Smooth, moulded acrylic will not pinch or bind tissue for greater patient comfort during examinations
- > The illumination system is completely sealed when inserted into the speculum
- > Easy one-handed operation provides multiple settings with positive locking positions

- 59005** KleenSpec 590 Series Disposable Vaginal Specula with Smoke Evacuator Tube, Small (case of 48)
- 59006** KleenSpec 590 Series Disposable Vaginal Specula with Smoke Evacuator Tube, Medium (case of 48)

