



Expectoration for adults and children

Probably the most popular mucus clearance devices used by professionals, patients and consumers with over 1 million sold worldwide.

Shaker Medic PLUS₊ for Hospitals / Clinics

Simple relief from excessive mucus and its associated problems

- Autoclavable to 134°C (273°F), multi-user, reuseable, easy to use and affordable
- Improved air flow for higher intensity vibrations, produces positive expiratory pressure (PEP) to open the airways
- Increases expectoration, improves lung function and reduces dyspnoea
- For conditions such as COPD, Cystic Fibrosis, bronchiectasis and bronchitis
- Two designs of user friendly mouthpieces included to assist older and younger patients



Shaker deluxe

Personal Mucus Clearance Device

- Respiratory device for mobilising pulmonary secretions such as mucus and catarrh
- Can be used for chronic conditions such as COPD, asthma, emphysema and acute problems like chesty coughs, flu and bronchitis
- Uses the latest innovative design technology
- Produces vibrations in the chest cavity to 'shake' stubborn mucus loose
- Aids expectoration, increases lung efficiency, helping to reduce breathlessness and fatigue
- Two designs of user friendly mouthpieces included to assist older and younger users
- Affordable, easy to use single user model which can be used at home or away from home



Shaker classic

Personal Mucus Clearance Device

- The original Shaker device offers simple and convenient relief from excessive mucus
- Can be used for chronic conditions such as COPD, asthma, emphysema and acute problems like chesty coughs, flu and bronchitis
- Affordable, easy to use single user model which can be used at home or away from home



For more information and specifications on the Shaker range from POWERbreathe visit powerbreathe.com