

Mobile stand for cable remote displays of seca scales and measuring rods seca 472



The stable solution for cable remote displays


This one is not so easy to bowl over: the stand seca 472 is the transportable solution for holding the display and control element of the digital scales and measuring rods seca 862, 884, 883, 242, 635, 634, 675, 674, 985 and 606.

Firm hold and easy to transport

The stand seca 472 is not only a practical alternative to the table support or wall fixture. Equipped with non-slip rubber feet, it also replaces the former model on transport castors and can be positioned everywhere quickly and securely. And its low weight of just three kilograms makes it easy to move around.

Sturdy, durable construction

The stable stand is made of high grade steel tubing and is thus of durable quality. The display and operating control element is held firmly by a robust clamp. The height and inclination angle of the display element can be adjusted to individual requirements. It can therefore be operated comfortably both from a seated or standing position.



Technical data

seca 472

Dimensions (WxHxD): 50 x 100 x 50 cm

Weight: 3 kg