



AsPEKT A703/712 digital recorders with LCD screen



AsPEKT A703/712 recorders:

- ▶ presentation of the ECG signal at the LCD screen,
- ▶ easily detachable patient cable,
- ▶ SD card support,
- ▶ Bluetooth connection with PC.

AsPEKT A703/712 digital recorders



AsPEKT A703/712 recorders:

- ▶ little portable dimensions,
- ▶ recording of 3 channels (AsPEKT 703)/ 12 channels (AsPEKT 712),
- ▶ SD card support,
- ▶ Bluetooth connection with PC.


Disclaimer: all materials, descriptions and images included in this brochure are for informative purposes only. Product images may differ from actual product appearance.

10.2017



HolCARD 24W

HolCARD 24W is an easy to operate software enabling efficient moving between ECG registration and analysis. In cooperation with AsPEKT recorders it creates a system to be installed on an user's PC and is dedicated to carry out both quick, basis and detailed, advanced analysis of ECG Holter examination.



- ▶ estimation of arrhythmias and of ST segment
- ▶ analyses of ECG registration in prospective and retrospective modes
- ▶ detection, reviewing and compiling quantitative reports within the scope of: VE, single-SVE, PSVT and many others
- ▶ results' presentation in a trend format, likewise ST trends in 3-D option
- ▶ analyses of arrhythmias of ventricular and supraventricular origins in registrations till 48 hours
- ▶ easy edition of all ventricular and supraventricular impulses owing to the used class editor function
- ▶ independent analyses of ST segment of each ECG channel for chosen by a user ST points of reference; J point as well as R-R base line

Disclaimer: all materials, descriptions and images included in this brochure are for informative purposes only. Product images may differ from actual product appearance.

10.2017

