

# Sun Autoclave N Class



## COMMODITIES PARAMETERS

Voltage	AC 220V/110v ±10%
Frequency	50/60Hz
Wire plug	GB 3 cores/EN 3 cores/ American 3 cores
Max Hp	2000VA
Tank Capacity	23L
Class	European Class N standard
Sterilizing temperature	134°C/121°C
Inside Dimension	Φ249×450mm
Net weight	39kg

## Product details

Model	Power	Frequency	Tank Capacity	Appearance size	Inside Dimension	Outside Dimension	G"W
SUN23-I (N)	2000W	50/60Hz±1Hz	23Litres	445X640X395mm	Ø249x450mm	740×550×500mm	39kg
SUN18-I (N)	1800W	50/60Hz±1Hz	18Litres	445X550X395mm	Ø249x355mm	670×550×500mm	33kg
SUN16-I (N)	1800W	50/60Hz±1Hz	16Litres	445X550X395mm	Ø230x360mm	670×550×500mm	32kg

## Technical Parameters

Sterilization level	European Class N standard
Sterilizing temperature	134°C/121°C
Display	Widescreen LED
Safety	Safety valve, Duble Pressure lock Automatic check.
Instrument containers	3 levels with 3 Round trays

## PRODUCT CHARACTERISTICS

1. Based on Europe Class N standard, installed with wrapped and unwrapped system, both of them have the choice on 121 and 134°C
2. Sterilizing were controlled by computer, most suitable for the dept of Stomatology, Ophthalmology, Surgery and Lab.
3. With safety system like: safety valve, manual door lock, error self-test and overheat protector.
4. The chamber were made by special Stainless Steel #304, working pressure from-0.9bar to 2.3 bar.
5. The bracket installed with 3 sterilizing plate.
6. Quick sterilizing process in 15-30 minutes(depend on the temperature)