

## UROFLOWMETER WIRELESS FLOWLINE CONNECT **NEW PRODUCTS**



Uroflussimetro with proprietary technology Wireless, range from 50-80 but closed doors, easy to use and portable dedicated in his briefcase.

### **Uroflowmeter Wireless Flowline Connect**

latest production and novelty of 2015. The Flowline Connect represents and meets the needs of all professionals who want to have a versatile space-saving and generation. Technology Touch Screen Tablet Portable 7.10.12 "high definition

Wireless sensor technology EV.SERVICE ITALIA owner, positioning up to 50-80 meters away. The transducer is compact ergonomic plastic can be cleaned easily and completely dismantled to be transported in its case dedicated. Operation even in battery to always have the tool in any situation.

The Flowline Connect is simple to use, the software functions are automated and make its use completely accessible even to those who are not computer experts.

The Flowline Connect has a transducer, which always keeps the calibration over time and easily controllable.

Stick connect proprietary, turn on the tablet and you're ready for your exam. The Report is customizable and can be printed on various media, from classical laser printer, inkjet up to those with portable thermal paper.

THE FLOWLINE CONNECT IS PROVIDED IN THE FOLLOWING CONFIGURATION:

TABLET TOUCH SCREEN 7"-10"

KEY WIRELESS TECHNOLOGY

SOFTWARE PRE-UPDATE Connect

PLASTIC FUNNEL PORT TRANSMITTER FROM 1 litre and BEKER (support for OPTIONAL HEIGHT ADJUSTMENT)

**12V POWER SUPPLY (for operation and BATTERY CHARGING)**

**Operator's Manual in English**

CASES for TRANSPORTING

EV.SERVICE ITALIA SRL

VIA UNGHERI 7/A 40012 CALDERARA DI RENO ( BOLOGNA ) TEL 051 728836 E-MAIL :Info@evserviceitalia.org

