

# EMERGENCY TROLLEYS

## EMERGENCY TROLLEYS - small and standard

Professional trolley with painted steel structure and external plastic finish equipped with 4 drawers, one lateral waste container, one extra lateral table, I.V. stand, defibrillator holder (plexiglass tray), oxygen holder up to 5 litres (back), 2 antistatic rubber castors.

- 27479 EMERGENCY TROLLEY - small 67x64xh 100 cm
- 27471 DIVIDER KIT for 27479 (drawer 45 cm)
- 27482 EMERGENCY TROLLEY - standard 82x64xh 100 cm
- 27472 DIVIDER KIT for 27482 (drawer 60 cm)

## DISPOSABLE SEALS

- 27473 DISPOSABLE SEALS - box of 100
- 27474 DISPOSABLE NUMBERED SEALS - box of 100

## EMERGENCY STANDARD



Transparent cardiac board

3 drawers  
60x48x13 cm

1 drawer  
60x48x26 cm

Built in perimetral protection



# KS MULTIFUNCTIONAL TROLLEYS

## MULTIFUNCTIONAL TROLLEYS: BASIC, DRESSING, EMERGENCY

Professional trolleys in 2 sizes (small and standard) and 3 models (Basic, Dressing and Emergency).

- antimicrobial epoxy painted body with seal lock (Emergency) or key lock (Basic and Dressing)
- top surface with 4-sided rim
- top surface with built-in compartment
- ball-bearing full-extension sliding drawer with self-closing system

- built-in handle
- 4 swivelling castors Ø 125 mm - 2 antistatic (only 45688/9) and with brake
- lateral waste container 15 l with lid and catheter holder (Dressing and Basic only)
- size: small 760x714xh 1,022 mm  
standard 911x704xh 1,075 mm

Lateral side with waste container 15 l with lid and catheter holder

45680



Dispenser with pole including 10 tilting bins

45685



Back with large oxygen tank compartment and stainless steel holding rail with hooks  
Lateral side with waste containers 15 l with lid

45689



Height adjustable I.V. pole with 2 hooks

Defibrillator/monitor rotating tray

45689



45681

2 small drawers: 600xh 75 mm  
2 medium drawers: 600xh 151 mm  
1 large drawer: 600xh 228 mm



Writing surface

Lateral container unit with 4 supply trays

45684

3 medium drawers: 450xh 151 mm  
1 large drawer: 450xh 228 mm



Lateral container unit with 4 supply trays

- 45680 BASIC TROLLEY - small
- 45681 BASIC TROLLEY - standard

- 45684 DRESSING TROLLEY - small
- 45685 DRESSING TROLLEY - standard  
Dressing trolley is supplied with upper dispenser (10 tilting bins).

- 45688 EMERGENCY TROLLEY - small
- 45689 EMERGENCY TROLLEY - standard  
Emergency trolley is supplied with rotating defibrillator/monitor tray, height adjustable I.V. pole with 2 hooks.