

## Classic Human Skull Model, 3 part - 3B Smart Anatomy

### Product Data



Weight	0.8 kg
Dimensions	20 x 13.5 x 15.5 cm
Brand	3B Scientific
Item No.	A20 [1020159]

- This product carries a 5 year warranty.

### Description

This 3 part medical quality human skull replica is a first choice for basic anatomical studies of the skull. It also makes a great present for medical professionals and students of medicine and allied health professions.

- High quality original human skull cast
- Handmade from hard, unbreakable plastic
- Highly accurate representation of the fissures, foramina, processes, sutures etc.
- Can be disassembled into skull cap, base of skull and mandible
- Mandible is mounted on a spring to easily demonstrate natural movement
- As an option, a 5 part brain can be purchased separately (C 18) and fits into the skull
- Each tooth individually cast for realistic dentition
- Now with magnetic connections

This quality skull model is a great addition to any study of the human skeletal system.

