

PRODUCT SPECIFICATION

REF 390104 FIRST AID KIT DIN



1 GENERAL PRODUCT DESCRIPTION

Description

A First Aid kit containing first aid products according to DIN 13157, is to be used in emergency situations (bleeding, breathing, burns, sprains, circulatory shock, minor eye injuries). The strong and durable case also has a Salvequick plaster dispenser on the outside allowing easy access to plasters without opening the case. Inside the well-organised case, all products are easily available and protected against dust, dirt, and moisture. The bright green and yellow coloured kit are clearly visible at a distance. Space inside for additional products.

The Medical device products inside are CE-labelled.

Accessory: Wall bracket, REF 5100008

Storage information

Dry storage at room temperature. For maximum protection, keep the products in their original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit (case)

<u>Product/Material</u>	<u>Qty</u>	<u>REF (refill)</u>	<u>Product/Material</u>	<u>Qty</u>	<u>REF (refill)</u>
First Aid Instructions	1	-	Eye Dressing Pad	2	-
First Aid instruction DIN	1	-	Salvequick Textile Plasters x 40	1	6444
Content List	1	-	Medical Face Mask	2	-
Cederroth Breathing Mask	1	2596	Salvequick Plastic Plasters x 45	1	6036
Gloves, pairs	2	2596	Salvequick Mixed Plasters	42	-
Cederroth Safety Hand Cleanser	4	2596	19x 72 mm	6	-
Cederroth Emergency Blanket	1	1892	25x72 mm	12	-
Scissors	1	270002	Finger bandage	6	-
Tweezers	1	1875	10x6 cm	12	-
Plastic Bag	2	-	Fingertip	6	-
Cederroth 4 in 1 Bloodstopper	3	1910	Plaster Tape Sensitive 5m x 2.5cm	1	310366
Cederroth 4 in 1 Bloodstopper mini	2	1911	Elastic bandage, 8 cm x 4 m	2	1882
Dressing 10 x 10 cm, sterile	6	-	Elastic bandage, 6 cm x 4 m	2	-
Dressing Cloth 60 x 80 cm	1	-	Salvequick Cold Pack	1	219600
Dressing Cloth 20 x 30 cm	5	-	Triangular Bandage	2	1889

Material

Case: Case: HDPE, handle and plaster dispenser: PP, latches: POM, button lock and inside dividers: PS, Inside cover: acrylic glass

First Aid products: Refer to product specifications for each product

Shelf life: 5 years from the production date for each component

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310613901049	7310803901040	7312023901046
Content:	1 kit	2 CU	6 SKU/layer x 4 layers
Size (w x l x h mm):	420 x 118 x 300	430 x 305 x 249	Pallet 800 x 1200 x 1146
Gross weight:	3409g	7150 g	197 kg
Packing material:	Carton 47 g	Corrugated board 332 g	-



Version: 05
Issued: 2022-05-19